


## DURABLE AND GENTLE ENTRAP POLYP RETRIEVAL NET

The Entrap™ retrieval system delivers a quick and efficient retrieval solution for polyps and tissue for histopathologic review. Utilizing advanced materials and production techniques, the Entrap facilitates improved consistency and control of the net deployment, position, and closure.



HANDLE MAY  
BE ROTATED

SMOOTH GLIDE

FIRM CATHETER



## KEY BENEFITS

### REFINED NET

Improved visibility through the net material for identification of the retrieval target. Durable enough for repeated deployment, but soft to gently capture and preserve the specimen.

### ADVANCED PROCESSES

Cinching, tearing, and binding of the net material is reduced by the minimal net selvage and bonded net attachment. This increases efficiency and minimizes the time required to reliably capture specimens.

### SMOOTH AND STRONG

A stiffer catheter material improves device insertion into the endoscope and positioning of the net. The additional columnar strength allows for gentle movement of a specimen for optimal capture.

## SPECIFICATIONS

### ENTRAP RETRIEVAL SYSTEM

Order Number	Description	Loop Diameter (mm)	Sheath Diameter (Fr)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
NT54511	Entrap Polyp Retrieval Net	30	7	230	2.8	10/Box
NT54515	Entrap Polyp Retrieval Net	30	7	230	2.8	2/Box

## OTHER PRODUCTS

### COLD SNARES EFFICIENT RESECTION



### BLUE BEACON TRUE BLUE GI MARKER



### POLYP TRAP SINGLE & MULTI-CHAMBER



### CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.



2855 Boardwalk Drive  
Ann Arbor, MI 48104  
877.552.4027

[MTEndoscopy.com](http://MTEndoscopy.com)  
[customer care@mtendoscopy.com](mailto:customer care@mtendoscopy.com)

MPW54000 Rev. 5