

RiteBite™

TWIN-WIRE BIOPSY FORCEPS

Twin-wire mechanism and precision stamping technology for assured tissue sampling.

- Pivoting jaws to aid in sampling where space may be limited
- Elongated jaw for greater tissue retention
- With-needle option for multiple bites in a single pass
- Entry - Exit indicator markings on sheath
- Flexible and pushable smooth catheter



KEY BENEFITS

TWIN-WIRE DESIGN

The twin-wire design allows the forcep head to rotate radially, while providing suitable closure force to the jaws, which may improve sampling ability and complete jaw closure.

ERGONOMIC HANDLE

Three-ring comfort with cushioned lining and short throw provides improved grip and accommodates smaller hands.

SUITABLE SAMPLE SIZE

The high yield family of forceps from Micro-Tech helps the physician to achieve a suitable amount of tissue for accurate and thorough diagnosis. The large capacity, jumbo and jumbo max forceps are designed to provide optimal bite along with spike options for multiple bites in a single pass.



SPECIFICATIONS

TWIN-WIRE BIOPSY FORCEPS STERILE

Order Number	Description	Jaw Style	Sheath Color*	Working Length (cm)	Minimum Channel Size (mm)	Package Units
BF40286	Large Capacity Forceps	Alligator	Blue	230	2.8	50/box
BF40306	Large Capacity Forceps with Spike	Alligator	Blue	230	2.8	50/box
BF40461	Jumbo Biopsy Forceps	Serrated	Blue	230	2.8	10/box
BF40451	Jumbo Biopsy Forceps with Spike	Serrated	Blue	230	3.2	10/box
BF40501	Jumbo Max Biopsy Forceps	Serrated	Blue	230	3.7	10/box
BF40491	Jumbo Max Biopsy Forceps with Spike	Serrated	Blue	230	3.7	10/box

*Subject to Change

CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.