

Sampling made simple

TWIN-WIRE BIOPSY FORCEPS

Biopsy forceps are the gold standard for tissue collection during endoscopic procedures, providing the accuracy, reliability, and consistency clinicians need. Micro-Tech's Twin-Wire Biopsy Forceps portfolio offers a comprehensive range of sizes to suit various procedural needs. Engineered with precision-stamped, twin-wire technology, these forceps deliver firm, consistent jaw closure, helping to achieve a reliable bite with every use.



KEY BENEFITS

TWIN-WIRE DESIGN

Our twin-wire mechanism is designed to support radial jaw rotation and firm closure, helping clinicians achieve accurate and consistent specimen collection, even in anatomically challenging cases.

ERGONOMIC HANDLE

Designed with clinician comfort in mind, our short-throw, cushioned three-ring handles was designed to improve grip helping to control, and minimize fatigue during high-volume procedures.

WIDE RANGE OF JAW STYLES

Whether you need standard, large, or jumbo capacity, with or without needles, our diverse lineup helps support optimal sampling based on lesion size, tissue type, and procedural goals.



RITEBITE™ BIOPSY FORCEPS STERILE

Order Number	Reference Number	Description	Jaw Style	Sheath Color	Working Length (cm)	Minimum Channel Size (mm)	Package Units
BF40661	SCNBF13-10022180	Standard Capacity Forceps	Alligator	Yellow	180 cm	2.8	10/Box
BF40641	SCNBF13-10122180	Standard Capacity Forceps with Spike	Alligator	Yellow	180 cm	2.8	10/Box
BF40671	SCNBF13-10122230	Standard Capacity Forceps	Alligator	Orange	230 cm	2.8	10/Box
BF40676	SCNBF13-10122230	Standard Capacity Forceps	Alligator	Orange	230 cm	2.8	50/Box
BF40651	SCNBF13-10122230	Standard Capacity Forceps with Spike	Alligator	Orange	230 cm	2.8	10/Box
BF40656	SCNBF13-10122230	Standard Capacity Forceps with Spike	Alligator	Orange	230 cm	2.8	50/Box
BF40281	LCNBF13-10024230	Large Capacity Forceps	Alligator	Orange	230 cm	2.8	10/Box
BF40286	LCNBF13-10024230	Large Capacity Forceps	Alligator	Orange	230 cm	2.8	50/Box
BF40301	LCNBF13-10124230	Large Capacity Forceps with Spike	Alligator	Orange	230 cm	2.8	10/Box
BF40306	LCNBF13-10124230	Large Capacity Forceps with Spike	Alligator	Orange	230 cm	2.8	50/Box
BF40461	JSNBF63-10028230	Jumbo Forceps	Serrated	Orange	230 cm	3.2	10/Box
BF40451	JSNBF63-10128230	Jumbo Forceps with Spike	Serrated	Orange	230 cm	3.2	10/Box
BF40491	JSNBF63-10133230	Jumbo Max Forceps with Spike	Serrated	Orange	230 cm	3.7	10/Box

CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.