

## DURACLIP<sup>™</sup> HEMOSTASIS CLIP

# DEPLOY WITH PRECISION



Now part of the Micro-Tech portfolio, DuraClip<sup>™</sup> continues to deliver the trusted performance and quality clinicians expect.

The DuraClip Hemostasis Clip delivers precise, reliable tissue approximation across a range of GI applications. Designed for control and simplicity, DuraClip can be opened and closed an unlimited number of times before deployment, helping clinicians reposition with confidence.

With responsive rotation for accurate positioning, flexible access in challenging anatomy, and intuitive deployment, DuraClip is designed to support efficient, controlled clip placement.

## KEY FEATURES

### CONTROLLED CLIP PLACEMENT

True 1:1 rotation with unlimited open and close capability before deployment supports precise positioning.

### ACCESS CHALLENGING ANATOMY

Flexible distal end enables rotation over ERCP scope elevator and accommodates insertion and rotation in retroflexed scope.

### DESIGNED FOR PROCEDURAL EFFICIENCY

Small stem profile and intuitive deployment support visualization and ease of use during clip placement.

## DURACLIP HEMOSTASIS CLIP

| UPN/Order Number | REF/Catalog Number | Opening Width (mm) | Sheath Diameter (mm) | Working Length (cm) | Minimum Channel Size (mm) | Package Units |
|------------------|--------------------|--------------------|----------------------|---------------------|---------------------------|---------------|
| RC40011          | DC0235             | 11                 | 2.6                  | 235                 | 2.8                       | 10/Box        |
| RC40031          | DC0235W            | 16                 | 2.6                  | 235                 | 2.8                       | 10/Box        |

### CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.