
MICRO-TECH® ENDOSCOPY

EXCEPTIONAL ENDOSCOPIC SOLUTIONS



MAXIMUM FLEXIBILITY

BILIARY STONE RETRIEVAL BALLOON

Micro-Tech's stone extraction balloon is comprised of a high-grade latex for reliability during endoscopic bile stone extraction. The 3 stage balloons come in above and below injection options to meet a wide range of clinical needs. The catheter features a 3-lumen design for insufflation, guidewire access, and contrast injection.

KEY BENEFITS

PROCEDURAL FLEXIBILITY
3-stage balloons are sized in 3 mm increments between 9 and 21 mm.

INSERTION AND TRACKING
The tapered atraumatic tip helps to ease insertion and improve tracking along the wire into biliary anatomy.

COMPATIBILITY
The balloon catheter is compatible with 0.035" guidewires.

STONE RETRIEVAL BALLOONS

STERILE EO

Order Number	Ref/ Catalog Number	Catheter Diameter (Fr)	Balloon Stages	Injection Port Location	Working Length (cm)	Package Units
SB75921	RSRB-P-9/12/15-20	7	9/12/15 mm	Above (distal to balloon)	200	1/Box
SB75931	RSRB-P-12/15/18-20	7	12/15/18 mm	Above (distal to balloon)	200	1/Box
SB75941	RSRB-P-15/18/21-20	7	15/18/21 mm	Above (distal to balloon)	200	1/Box
SB76131	RSRB-P-9/12/15-20-B	7	9/12/15 mm	Below (proximal to balloon)	200	1/Box*
SB76141	RSRB-P-12/15/18-20-B	7	12/15/18 mm	Below (proximal to balloon)	200	1/Box*
SB76151	RSRB-P-15/18/21-20-B	7	15/18/21 mm	Below (proximal to balloon)	200	1/Box*

*Non-Stock Item

Rx Only. See Instructions for Use for compatible indications, contraindications, warnings, precautions, and instructions for use.

US.MTMED.COM
customercare@mtmed.com

2855 Boardwalk Drive
Ann Arbor, MI 48104
877.552.4027