MICRO-TECH®
ENDOSCOPY

EXCEPTIONAL ENDOSCOPIC SOLUTIONS

PRECISE CANNULATION

SPHINCTEROTOME

The Micro-Tech sphincterotome features a short and long wire mechanism in one device. The rotation wheel provides maximum control during cannulation and sphincterotomy.



KEY BENEFITS

FLEXIBILITY
Compatible with standard guidewires, the rotatable tip can be positioned to ensure proper cannulation and sphincterotomy positioning.

CONTROL
The c-channel catheter allows for a splitting movement to separate the wire from the catheter and the rotation wheel assists in fine tuned tip adjustments for successful physician guided cannulation.

ATRAUMATIC TIP
The tapered atraumatic tip helps to ease insertion and tracking.

SHOW PRELOADED SPHINCTEROTOMES STERILE EO

Order Number	Ref/ Catalog Number	Tip Length (mm)	Distal Tip O.D. (Fr)	Cutting Wire (mm)	Cutting Wire Insulation	Preloaded Wire Length (cm)	Preloaded Wire O.D. (in)	Working Length (cm)	Package Units
SS75761	RDSP-0520-221-1	5	3.9	20	Insulated	260	0.025	200	1/Box
SS75781	RDSP-0520-231-1	5	3.9	20	Insulated	450	0.025	200	1/Box*
SS16121	RDSP-0525-221-1	5	3.9	25	Insulated	260	0.025	200	1/Box
SS16141	RDSP-0525-231-1	5	3.9	25	Insulated	450	0.025	200	1/Box*
SS16181	RDSP-0530-221-1	5	3.9	30	Insulated	260	0.025	200	1/Box
SS16201	RDSP-0530-231-1	5	3.9	30	Insulated	450	0.025	200	1/Box*
SS75041	RDSP-0520-221	5	4.5	20	Insulated	260	0.035	200	1/Box
SS75061	RDSP-0520-231	5	4.5	20	Insulated	450	0.035	200	1/Box*
SS75101	RDSP-0525-221	5	4.5	25	Insulated	260	0.035	200	1/Box
SS75121	RDSP-0525-231	5	4.5	25	Insulated	450	0.035	200	1/Box*
SS75161	RDSP-0530-221	5	4.5	30	Insulated	260	0.035	200	1/Box
SS75181	RDSP-0530-231	5	4.5	30	Insulated	450	0.035	200	1/Box*

Preloaded with Hydro-Slide Straight Tip Guidewire
*Non-Stock Item

SHOW SPHINCTEROTOMES STERILE EO

Order Number	Ref/ Catalog Number	Tip Length (mm)	Distal Tip O.D. (Fr)	Cutting Wire (mm)	Cutting Wire Insulation	Recommended Wire O.D. (in)	Working Length (cm)	Package Units
SS75741	RDSP-0520-211-1	5	3.9	20	Insulated	0.025	200	1/Box
SS75801	RDSP-0525-211-1	5	3.9	25	Insulated	0.025	200	1/Box
SS16161	RDSP-0530-211-1	5	3.9	30	Insulated	0.025	200	1/Box
SS75021	RDSP-0520-211	5	4.5	20	Insulated	0.035	200	1/Box
SS75081	RDSP-0525-211	5	4.5	25	Insulated	0.035	200	1/Box
SS75141	RDSP-0530-211	5	4.5	30	Insulated	0.035	200	1/Box

Rx Only. See Instructions for Use for compatible indications, contraindications, warnings, precautions, and instructions for use.

US.MTMED.COM
customercare@mtmed.com
2855 Boardwalk Drive
Ann Arbor, MI 48104
877.552.4027