
MICRO-TECH®
ENDOSCOPY

EXCEPTIONAL ENDOSCOPIC SOLUTIONS

ENDOSCOPY
PRODUCT CATALOG

Partnering with physicians to deliver innovative, high-quality devices




EXCEPTIONAL **EXPERIENCE**

COST **CONTAINMENT**

PROCEDURAL **EFFICIENCY**

VALUE **INNOVATION**

At Micro-Tech Endoscopy, quality is the cornerstone of everything we do. We strive to deliver an exceptional experience by combining the highest quality state-of-the-art products with the support you can rely on to provide your patients with the finest care.



We partner with physicians and industry leaders to discover new ways to increase procedural efficiency. Our team works to develop and refine devices and services that aim to reduce procedure time, complications, or the need for additional therapies or treatments.

Here at Micro-Tech, cost containment means operating differently. We employ world-leading production and logistics technologies, automation, design, and lean practices to take cost out of our operations. We focus our resources on what matters, such as product and technological advancements and our customer's experience. As a result, we have a proven record of delivering value.

OUR PROMISE OF AN EXCEPTIONAL EXPERIENCE,
improved procedural efficiency, and cost containment, all combine to deliver our goal of value innovation. It is our focus and what drives us to innovate and discover more ways to create value that can make a meaningful difference in healthcare delivery and economics for everyone.

Micro-Tech Endoscopy's mission is to create quality, cutting-edge products, and services that help improve patient outcomes. We are committed to providing an exceptional experience and keeping costs low.

TABLE OF CONTENTS

BILIARY 4

BIOPSY 9

CHOLANGIOSCOPY 16

DILATION 23

EUS NEEDLES 29

HEMOSTASIS 33

INFECTION PREVENTION 43

NITINOL STENTS 55

EMR-ESD 60

POLYPECTOMY- RETRIEVAL 73

PRODUCT REFERENCE CHART 82

TERMS AND POLICIES 87

HOW TO ORDER



CALL

877.552.4027
Monday – Friday
8:30 AM – 6:00 PM EST



EMAIL

customercare@mtmed.com



FAX

734.272.4160



MAIL

Micro-Tech Endoscopy USA
2855 Boardwalk Drive
Ann Arbor, MI 48104



INTERESTED IN PLACING ORDERS VIA EDI?

edi.mtam@mtmed.com

BILIARY

GUIDEWIRES

PRELOADED SPHINTEROTOMES

SPHINTEROTOMES

STONE RETRIEVAL BALLOONS



ACCELERATED ACCESS

HYDRO-SLIDE GUIDEWIRE

A reliable guidewire is the core of a successful procedure. The Hydro-Slide Guidewire provides the stiffness, visibility, and atraumatic tip needed to bridge strictures in the GI tract.

The PTFE jacket, Hydrophilic coating, and Nitinol core provide improved pushability, while the soft tungsten-loaded tip allows for high visibility under fluoroscopy. The yellow and blue striped jacket is an easy-to-see indicator when advancing the guidewire under direct visualization.

- Straight and angled tips
- 260 cm and 450 cm working lengths
- Individually packaged, 2 per box

HYDRO-SLIDE .025" GUIDEWIRE STERILE ② X-RAY

UPN/Order Number	REF/Catalog Number	Description	Working OD (in)	Working Length (cm)	Package Units
GW10011	MTN-BM-63/26-A	Hydro-Slide Guidewire Straight	.025	260	2/Box
GW10021	MTN-BM-63/26-A-J	Hydro-Slide Guidewire Angled	.025	260	2/Box
GW10031	MTN-BM-63/45-A	Hydro-Slide Guidewire Straight	.025	450	2/Box*
GW10041	MTN-BM-63/45-A-J	Hydro-Slide Guidewire Angled	.025	450	2/Box*


HYDRO-SLIDE .035" GUIDEWIRE STERILE ② X-RAY

UPN/Order Number	REF/Catalog Number	Description	Working OD (in)	Working Length (cm)	Package Units
GW10051	MTN-BM-89/26-A	Hydro-Slide Guidewire Straight	.035	260	2/Box
GW10071	MTN-BM-89/26-A-J	Hydro-Slide Guidewire Angled	.035	260	2/Box
GW10091	MTN-BM-89/45-A	Hydro-Slide Guidewire Straight	.035	450	2/Box*
GW10111	MTN-BM-89/45-A-J	Hydro-Slide Guidewire Angled	.035	450	2/Box*

HYDRO-SLIDE .035" STIFF GUIDEWIRE STERILE ② X-RAY

UPN/Order Number	REF/Catalog Number	Description	Working OD (in)	Working Length (cm)	Package Units
GW10061	MTN-BM-89/26-A-S	Hydro-Slide Guidewire Stiff Straight	.035	260	2/Box
GW10081	MTN-BM-89/26-A-J-S	Hydro-Slide Guidewire Stiff Angled	.035	260	2/Box
GW10101	MTN-BM-89/45-A-S	Hydro-Slide Guidewire Stiff Straight	.035	450	2/Box*
GW10121	MTN-BM-89/45-A-J-S	Hydro-Slide Guidewire Stiff Angled	.035	450	2/Box*

*Non-Stock Item



PRECISE CANNULATION

SPHINCTEROTOME

The Micro-Tech sphincterotome features a patent-pending[1] opening mechanism to facilitate both short-wire and long-wire procedures with one device. The rotatable design and locking mechanism provide maximum control during cannulation and sphincterotomy.

KEY BENEFITS

FLEXIBILITY
Compatible with 0.035" guidewires, the rotatable tip can be actuated to ensure proper cannulation and sphincterotomy positioning.

CONTROL
The c-channel catheter allows for a splitting movement to separate wire from catheter for physician-guided cannulation.

ATRAUMATIC TIP
The tapered atraumatic tip helps to ease insertion and tracking.


SHOW PRELOADED SPHINCTEROTOMES STERILE EO

Order Number	Ref/ Catalog Number	Tip Length (mm)	Distal Tip O.D. (Fr)	Cutting Wire (mm)	Cutting Wire Insulation	Preloaded Wire Length (cm)	Preloaded Wire O.D. (in)	Working Length (cm)	Package Units
SS75761	RDSP-0520-221-1	5	3.9	20	Insulated	260	0.025	200	1/Box
SS75781	RDSP-0520-231-1	5	3.9	20	Insulated	450	0.025	200	1/Box*
SS16121	RDSP-0525-221-1	5	3.9	25	Insulated	260	0.025	200	1/Box
SS16141	RDSP-0525-231-1	5	3.9	25	Insulated	450	0.025	200	1/Box*
SS16181	RDSP-0530-221-1	5	3.9	30	Insulated	260	0.025	200	1/Box
SS16201	RDSP-0530-231-1	5	3.9	30	Insulated	450	0.025	200	1/Box*
SS75041	RDSP-0520-221	5	4.5	20	Insulated	260	0.035	200	1/Box
SS75061	RDSP-0520-231	5	4.5	20	Insulated	450	0.035	200	1/Box*
SS75101	RDSP-0525-221	5	4.5	25	Insulated	260	0.035	200	1/Box
SS75121	RDSP-0525-231	5	4.5	25	Insulated	450	0.035	200	1/Box*
SS75161	RDSP-0530-221	5	4.5	30	Insulated	260	0.035	200	1/Box
SS75181	RDSP-0530-231	5	4.5	30	Insulated	450	0.035	200	1/Box*

Preloaded with Hydro-Slide Straight Tip Guidewire
*Non-Stock Item

SHOW SPHINCTEROTOMES STERILE EO

Order Number	Ref/ Catalog Number	Tip Length (mm)	Distal Tip O.D. (Fr)	Cutting Wire (mm)	Cutting Wire Insulation	Recommended Wire O.D. (in)	Working Length (cm)	Package Units
SS75741	RDSP-0520-211-1	5	3.9	20	Insulated	0.025	200	1/Box
SS75801	RDSP-0525-211-1	5	3.9	25	Insulated	0.025	200	1/Box
SS16161	RDSP-0530-211-1	5	3.9	30	Insulated	0.025	200	1/Box
SS75021	RDSP-0520-211	5	4.5	20	Insulated	0.035	200	1/Box
SS75081	RDSP-0525-211	5	4.5	25	Insulated	0.035	200	1/Box
SS75141	RDSP-0530-211	5	4.5	30	Insulated	0.035	200	1/Box



MAXIMUM FLEXIBILITY

BILIARY STONE RETRIEVAL BALLOON

Micro-Tech's stone extraction balloon is comprised of a high-grade latex for reliability during endoscopic bile stone extraction. The 3-stage balloons come in above and below injection options to meet a wide range of clinical needs. The catheter features a 3-lumen design for insufflation, guidewire access, and contrast injection.

KEY BENEFITS

PROCEDURAL FLEXIBILITY
3-stage balloons are sized in 3 mm increments between 9 and 21 mm.

INSERTION AND TRACKING
The tapered atraumatic tip helps to ease insertion and improve tracking along the wire into biliary anatomy.

COMPATIBILITY
The balloon catheter is compatible with 0.035 guidewires.

SHOW STONE RETRIEVAL BALLOONS STERILE EO

UPN/Order Number	REF/Catalog Number	Catheter Diameter (Fr)	Balloon Stages	Injection Port Location	Working Length (cm)	Package Units
SB75921	RSRB-P-9/12/15-20	7	9/12/15 mm	Above (distal to balloon)	200	1/Box
SB75931	RSRB-P-12/15/18-20	7	12/15/18 mm	Above (distal to balloon)	200	1/Box
SB75941	RSRB-P-15/18/21-20	7	15/18/21 mm	Above (distal to balloon)	200	1/Box
SB76131	RSRB-P-9/12/15-20-B	7	9/12/15 mm	Below (proximal to balloon)	200	1/Box*
SB76141	RSRB-P-12/15/18-20-B	7	12/15/18 mm	Below (proximal to balloon)	200	1/Box*
SB76151	RSRB-P-15/18/21-20-B	7	15/18/21 mm	Below (proximal to balloon)	200	1/Box*

*Non-Stock Item



SHOW:

GUIDE WIRE LOCKING DEVICE

The SHOW Short Wire Locking Device simplifies and enhances the process of guidewire management during ERCP procedures. By fixing the guidewire in place, the SHOW Short Wire Locking Device allows for seamless device exchanges while minimizing the risk of wire displacement.

GUIDE WIRE LOCKING DEVICE

UPN/Order Number	REF/Catalog Number	Description	Package Units
GL62781	MT-RGL-O	SHOW Guidewire Locking Device, Oly, Fuj	10/Box
GL62791	MT-RGL-P	SHOW Guidewire Locking Device, Pentax	10/Box

BIOPSY

RITEBITE BIOPSY FORCEPS

TRUBITE BIOPSY FORCEPS

GI CYTOLOGY BRUSHES

MULTICROC


PRECISOR EXL

Sampling made simple

RITEBITE™ BIOPSY FORCEPS



Biopsy forceps are the gold standard for tissue collection during endoscopic procedures, providing the accuracy, reliability, and consistency clinicians need. Micro-Tech's Twin-Wire Biopsy Forceps portfolio offers a comprehensive range of sizes to suit various procedural needs. Engineered with precision-stamped, twin-wire technology, these forceps deliver firm, consistent jaw closure, helping to achieve a reliable bite with every use.



KEY FEATURES

TWIN-WIRE DESIGN

Our twin-wire mechanism is designed to support radial jaw rotation and firm closure, helping clinicians achieve accurate and consistent specimen collection, even in anatomically challenging cases.

ERGONOMIC HANDLE

Designed with clinician comfort in mind, our short-throw, cushioned three-ring handles were designed to improve grip helping to control, and minimize fatigue during high-volume procedures.

WIDE RANGE OF JAW STYLES


Whether you need standard, large, or jumbo capacity, with or without needles, our diverse lineup helps support optimal sampling based on lesion size, tissue type, and procedural goals.

RITEBITE™ BIOPSY FORCEPS STERILE

Order Number	Reference Number	Description	Jaw Style	Sheath Color	Working Length (cm)	Minimum Channel Size (mm)	Package Units
BF40661	SCNBF13-10022180	Standard Capacity Forceps	Alligator	Yellow	180	2.8	10/Box
BF40641	SCNBF13-10122180	Standard Capacity Forceps with Spike	Alligator	Yellow	180	2.8	10/Box
BF40671	SCNBF13-10022230	Standard Capacity Forceps	Alligator	Orange	230	2.8	10/Box
BF40676	SCNBF13-10022230	Standard Capacity Forceps	Alligator	Orange	230	2.8	50/Box
BF40651	SCNBF13-10122230	Standard Capacity Forceps with Spike	Alligator	Orange	230	2.8	10/Box
BF40656	SCNBF13-10122230	Standard Capacity Forceps with Spike	Alligator	Orange	230	2.8	50/Box
BF40271	LCNBF13-10024180	Large Capacity Forceps	Alligator	Yellow	180	2.8	10/Box
BF40291	LCNBF13-10124180	Large Capacity Forceps with Spike	Alligator	Yellow	180	2.8	10/Box
BF40281	LCNBF13-10024230	Large Capacity Forceps	Alligator	Orange	230	2.8	10/Box
BF40286	LCNBF13-10024230	Large Capacity Forceps	Alligator	Orange	230	2.8	50/Box
BF40301	LCNBF13-10124230	Large Capacity Forceps with Spike	Alligator	Orange	230	2.8	10/Box
BF40306	LCNBF13-10124230	Large Capacity Forceps with Spike	Alligator	Orange	230	2.8	50/Box
BF40461	JSNBF63-10028230	Jumbo Forceps	Serrated	Orange	230	3.2	10/Box
BF40451	JSNBF63-10128230	Jumbo Forceps with Spike	Serrated	Orange	230	3.2	10/Box
BF40491	JSNBF63-10133230	Jumbo Max Forceps with Spike	Serrated	Orange	230	3.7	10/Box

BIOPSY

TRUBITE BIOPSY FORCEPS



With our wide selection of cost-effective disposable biopsy forceps, you are perfectly equipped for every procedure.

Available jaw diameters of 1.8 mm, and 2.3 mm are offered in Oval, Serrated, Alligator jaws, with or without a spike.

- Available lengths of 160 cm and 230 cm
- All forceps feature entry and exit indicator markings



BIOPSY FORCEP FOR 2.0 MM SCOPE CHANNELS STERILE

UPN/Order Number	REF/Catalog Number	Description	Jaw Style	Sheath Color	Working Length (cm)	Jaw O.D. (mm)	Package Units
BF40021	NBF02-11018160	Mini Forceps	Oval	Yellow	160	1.8	10/Box
BF40051	NBF02-11118160	Mini Forceps with Spike	Oval	Yellow	160	1.8	10/Box
BF40141	NBF12-11018160	Mini Forceps	Alligator	Yellow	160	1.8	10/Box
BF40401	NBF12-11118160	Mini Forceps with Spike	Alligator	Yellow	160	1.8	10/Box

BIOPSY FORCEP FOR 2.8 MM SCOPE CHANNELS STERILE

UPN/Order Number	REF/Catalog Number	Description	Jaw Style	Sheath Color	Working Length (cm)	Jaw O.D. (mm)	Package Units
BF40091	NBF03-11023230	Standard Biopsy Forceps	Oval	Blue	230	2.3	10/Box
BF40121	NBF03-11123230	Standard Biopsy Forceps with Spike	Oval	Blue	230	2.3	10/Box
BF40096	NBF03-11023230	Standard Biopsy Forceps	Oval	Blue	230	2.3	50/Box
BF40126	NBF03-11123230	Standard Biopsy Forceps with Spike	Oval	Blue	230	2.3	50/Box
BF40201	NBF13-11023230	Standard Biopsy Forceps	Alligator	Blue	230	2.3	10/Box
BF40231	NBF13-11123230	Standard Biopsy Forceps with Spike	Alligator	Blue	230	2.3	10/Box
BF40206	NBF13-11023230	Standard Biopsy Forceps	Alligator	Blue	230	2.3	50/Box
BF40236	NBF13-11123230	Standard Biopsy Forceps with Spike	Alligator	Blue	230	2.3	50/Box
BF40441	NBF63-11023230	Standard Biopsy Forceps	Serrated	Blue	230	2.3	10/Box
BF40421	NBF63-11123230	Standard Biopsy Forceps with Spike	Serrated	Blue	230	2.3	10/Box
BF40446	NBF63-11023230	Standard Biopsy Forceps	Serrated	Blue	230	2.3	50/Box
BF40426	NBF63-11123230	Standard Biopsy Forceps with Spike	Serrated	Blue	230	2.3	50/Box

BIOPSY FORCEP HOT STERILE

UPN/Order Number	REF/Catalog Number	Description	Jaw Style	Sheath Color	Working Length (cm)	Jaw O.D. (mm)	Package Units
HF48521	HBF55-11023230	Hot Biopsy Forceps Lower	Oval	Red	230	2.3	10/Box
HF48541	HBF65-11023230	Hot Biopsy Forceps Lower	Serrated	Red	230	2.3	10/Box



GI CYTOLOGY BRUSHES

Stiff brushes for reliable tissue sampling

CYTOLOGY BRUSH STERILE

UPN/Order Number	REF/Catalog Number	Description	Brush O.D. (mm)	Catheter O.D. (mm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
CY49021	CYB-20S-18C	Cytology Brush	2.0	2.3	180	2.8	20/Box
CY49051	CYB-30S-24C	Cytology Brush	3.0	2.3	240	2.8	20/Box

MULTI-SAMPLING BIOPSY ADVANTAGE

MULTICROC



Tissue collection during endoscopic procedures requires efficiency, reliability, and precision. Biopsies often require multiple passes to collect adequate tissue samples, leading to longer procedure times, higher resource use and variable sample quality.

MultiCROC® Single Use Multiple Sampling Biopsy Forceps are designed to eliminate these inefficiencies. With the ability to collect five high-quality biopsy samples in a single pass, MultiCROC® improves workflow, optimizes resource use, and supports consistent sampling for both upper and lower GI procedures.

KEY FEATURES

FIVE SAMPLES IN A SINGLE PASS
Collect five targeted biopsy specimens in one pass through the endoscope, minimizing insertion time while maximizing sample yield.

ALLIGATOR-STYLE JAWS – 2.4 MM DIAMETER
Precision-engineered jaws ensure excellent tissue grip and clean cuts, optimized for a minimum working channel of 2.8 mm.

PATENTED DESIGN WITH LARGE BITING FORCE
Delivers superior tissue quality with defined borders and minimal crush artifact, enhancing the consistency of the samples.

OPTIMIZED JAW GEOMETRY
Facilitates efficient specimen release into fixative, streamlining tissue handling and preserving sample integrity for laboratory analysis.

MULTICROC STERILE

UPN/Order Number	REF/Catalog Number	Description	Package Units
MC924A-10	MC924A-10	MultiCROC® Biopsy Forceps, 240 cm, Multibite, 2.4 mm O.D.	10/Box

PRECISOR® EXL™
DISPOSABLE BIOPSY FORCEPS

CUT WITH CONFIDENCE



Now part of the Micro-Tech portfolio, Precisor® EXL™ continues to deliver the trusted performance and quality clinicians expect.

Precisor EXL Disposable Biopsy Forceps are designed to cut tissue, not tear, helping clinicians obtain high-quality specimens with confidence. Featuring CAM jaw technology that eliminates exposed linkage arms and hinges, Precisor EXL supports precise tissue acquisition while promoting smooth passage through the endoscope channel.

PRECISOR EXL BIOPSY FORCEPS

UPN/Order Number	REF/Catalog Number	Description	Jaw O.D (mm)	Jaw Style	Sheath Color	Working Length (cm)	Minimum Channel Size (mm)	Package Units
BF51811	000381	Standard Capacity Biopsy Forceps	2.3	Alligator	Yellow	160	2.8	10
BF51821	000382	Standard Capacity Biopsy Forceps with Spike	2.3	Alligator	Yellow	160	2.8	10
BF51831	000383	Standard Capacity Biopsy Forceps	2.3	Oval	Yellow	160	2.8	10
BF51841	000384	Standard Capacity Biopsy Forceps with Spike	2.3	Oval	Yellow	160	2.8	10
BF51851	000385	Standard Capacity Biopsy Forceps	2.3	Alligator	Blue	230	2.8	10
BF51861	000386	Standard Capacity Biopsy Forceps with Spike	2.3	Alligator	Blue	230	2.8	10
BF51881	000389	Standard Capacity Biopsy Forceps	2.3	Oval	Blue	230	2.8	10
BF51871	000388	Standard Capacity Biopsy Forceps with Spike	2.3	Oval	Blue	230	2.8	10
BF63851	00389B	Standard Capacity Biopsy Forceps	2.3	Oval	Blue	230	2.8	50
BF63831	00386B	Standard Capacity Biopsy Forceps with Spike	2.3	Alligator	Blue	230	2.8	50
BF52071	00821C	Jumbo Biopsy Forceps with Spike	3.1	Oval	Blue	230	3.2	10
BF52091	00823C	Jumbo Biopsy Forceps with Spike	3.1	Alligator	Blue	230	3.2	10
BF52141	100512	Pediatric Biopsy Forceps	1.8	Oval	Blue	160	2	10
BF52151	100513	Pediatric Biopsy Forceps	1.8	Alligator	Blue	160	2	10

CHOLANGIOSCOPY

EYEMAX SINGLE-USE CHOLANGIOSCOPE

EYEMAX EXTRACTION BASKET


EYEMAX 650 SERIES SCOPES

EYEMAX SLIM

WALZ EHL

NEXT-GEN VISUALIZATION

EYEMAX SINGLE-USE CHOLANGIOSCOPE



Discover the advantages of enhanced resolution by using the eyeMAX™ Cholangioscope. With superior optical clarity, 120° field of view and improved flexibility, eyeMAX provides unmatched access and clear visualization to aid diagnosis. Available in two sizes, the 3.9 mm O.D. scope includes a 2.0 mm working channel that accommodates our 1.6 mm eyeMAX biopsy forceps, potentially increasing biopsy yield. Meanwhile, the 3.2 mm version comes equipped with a 1.2 mm working channel that works seamlessly with the 1.0 mm eyeMAX biopsy forceps, enabling deeper access into the patient's anatomy.



eyeMAX™

CHOLANGIOSCOPES STERILE

UPN/Order Number	REF/Catalog Number	Version	Insertion Diameter (mm)	Length (cm)	Working Channel Diameter (mm)	Package Units
CD63681	CDS11007	H - Olympus only	3.9	220	2.0	1/Box
CD63701	CDS11008	C - Universal fitment	3.9	220	2.0	1/Box
CD63691	CDS22007	H - Olympus only	3.2	220	1.2	1/Box
CD63761	CDS22008	C - Universal fitment	3.2	220	1.2	1/Box

DIGITAL CONTROLLER

UPN/Order Number	REF/Catalog Number	Voltage	Frequency	Functions
SW71212	BS-W-100	120V	60Hz	LED, white color balance, brightness, zoom, image shape change


FORCEPS STERILE

UPN/Order Number	REF/Catalog Number	Description	Length (cm)	O.D. (mm)	Package Units
DF49425	BF10006	Biopsy Forceps, eyeMAX, Cam Alligator	290cm PTFE	1	2/Box
DF49365	BF16006	Biopsy Forceps, eyeMAX, Twin Wire Alligator	290cm PE	1.6	2/Box

CLEAR VISION
CONFIDENT STONE REMOVAL

EYEMAX EXTRACTION BASKET

Designed for use with our family of eyeMAX™ cholangioscopes, the eyeMAX extraction basket features a rotatable four-wire spiral design allowing for optimal positioning to effectively and efficiently capture and remove stones, stone fragments, and foreign bodies while under direct visualization.



KEY BENEFITS

FOUR WIRE SPIRAL DESIGN

Enhances gripping strength and flexibility, allowing for secure capture and easy manipulation.

REDUCED RADIATION EXPOSURE

By performing procedures under direct visualization, the eyeMAX extraction basket reduces the need for fluoroscopy thus minimizing radiation exposure.

PRECISE ROTATION

The 4-wire spiral design of the extraction basket, combined with enhanced rotation control, allows for precise maneuverability and secure capture and release of stones and foreign bodies.

EYEMAX™ EXTRACTION BASKET


UPN/Order Number	REF/Catalog Number	Description	Basket Diameter (mm)	Working Length (cm)
BE11721	CEB00000	eyeMAX Extraction Basket 4 Spiral 10 mm x 30 mm	10 mm	290 cm
BE11731	CEB01010	eyeMAX Extraction Basket 4 Spiral 15 mm x 40 mm	15 mm	290 cm




NEXT-GEN VISUALIZATION

EYEMAX 650 SERIES SCOPES

Designed for minimally invasive surgical biliary procedures, the 650 series eyeMAX scope provides practitioners with the tools they need for improved efficiency, visualization, and control. With crystal-clear imaging, eyeMAX allows for better diagnostic accuracy and facilitates complex therapeutic interventions. The 650 series eyeMAX scope is offered in two scope sizes: 3.9 mm OD and a 3.2 mm OD to provide the flexibility needed for comprehensive common bile duct exploration and treatment.





SCOPES

UPN/Order Number	REF/Catalog Number	Insertion Diameter (mm)	Length (cm)	Working Channel Diameter (mm)	Package Units
CD63241	CDS11005	3.9	65	2.0	1/Box
CD63251	CDS22005	3.2	65	1.2	1/Box

STONE RETRIEVAL BASKET

UPN/Order Number	REF/Catalog Number	Description	Basket Diameter (mm)	Working Length (cm)
BE11861	CEB00003	eyeMAX Extraction Basket 4 Spiral 10 mm x 30 mm	10	120
BE11871	CEB01013	eyeMAX Extraction Basket 4 Spiral 15 mm x 40 mm	15	120


NAVIGATE NARROWER. REACH FURTHER.
CHOLANGIOSCOPY REDEFINED.


EYEMAX SLIM


The next-gen eyeMAX™ Slim combines a streamlined, ultra-slim design with high-quality visualization to support direct biliary access in a wide range of cholangioscopy procedures. Designed for improved maneuverability through complex anatomy and compatibility with commonly used therapeutic duodenoscopes, eyeMAX Slim empowers physicians with the control, access, and image clarity needed to confidently diagnose and treat biliary disease.



eyeMAX™ Slim







With high-resolution imaging and compatibility with commonly used EHL systems and accessories, eyeMAX Slim supports precise diagnosis and effective treatment throughout the biliary and pancreatic tracts, advancing the capabilities of cholangioscopy.

KEY FEATURES



6.5 FR
0.75 MM WC



7.5 FR
1.2 MM WC



9.8 FR
2.0 MM WC




EYEMAX SLIM

UPN/Order Number	REF/Catalog Number	Description	Insertion Diameter (Fr)	Maximum O.D. (mm)	Working Channel Diameter (mm)	Package Units
AM71601	CDS33101	C- Universal Fitment	6.5	2.4	0.75	1/Box
AM71611	CDS44101	C- Universal Fitment	7.5	2.8	1.2	1/Box
AM71631	CDS66101	C- Universal Fitment	9.8	3.6	2.0	1/Box

SIMPLE & EFFECTIVE FOR LARGE STONE MANAGEMENT

EHL LITHOTRON EL27 AND ACCESSORIES

The Walz Lithotron EL27 Compact offers a secure, easy, and convenient way to administer bipolar ElectroHydraulic Lithotripsy (EHL).



KEY BENEFITS

CONTROL
Offering 3 different intensity options (up to 900 mJ), allowing users the ability to adjust energy levels based on stone densities.

SECURITY
Reduced tensile phases of concentrated energy enable accurate stone fragmentation while minimizing harm to soft tissues.

EFFICIENCY
The centered micro-shockwaves delivered by the probe facilitates targeting with fast pulses and sharply rising amplitudes.

EHL SYSTEM

UPN/Order Number	REF/Catalog Number	Description	Package Units
EL27C115VEHL	EL27C115VEHL	Lithotriptor EL-27 Compact 115V EHL main unit incl. Probe cable, foot switch, power cord.	1/box

EHL PROBES

UPN/Order Number	REF/Catalog Number	Probe Diameter (Fr)	Probe Diameter (mm)	Working Length (cm)	Package Units
P243000FST	P243000FST	2.4	0.8	300	1/box
P303000FST	P303000FST	3.0	1.0	300	1/box
P453000FST	P453000FST	4.5	1.5	300	

*Note 3.0 Fr probe can be used with all 3.2 mm & 3.9 mm scope but recommended for 3.9 mm eyeMAX scope with 2.0 mm channel.

**Note 4.5 Fr probe to be used exclusively with 3.9 mm eyeMAX scope with 2.0 mm channel

DILATION

MULTI-STAGE DILATION BALLOONS

ABC 3-STEP DILATION BALLOONS


INFLATION DEVICE

DILATION ACCESSORIES

STRENGTH AND CONTROL

MULTI-STAGE DILATION BALLOONS

Choose strength and control when dilating strictures with a dilation balloon.



Reliant® and Rapide® Multi-Stage balloons feature a premium clear Pebax® balloon with rounded shoulders for balloon endoscopy.

Multi-stage balloons allow for inflation to 3 distinct sizes to allow multi-step dilation in one balloon, providing controlled dilation and eliminating the need for multiple balloons in a procedure.

Fixed wire esophageal balloons offer fast, reliable performance.

Preloaded .035 wire-guided balloons assist placement in tortuous or altered anatomy.



.035 COMPATIBLE

KEY BENEFITS

CHOICE OF MATERIAL

The high-quality Pebax material produces a balloon with the features and performance your providers demand. Accurate three-size radial expansion, rounded shoulders, clear visibility and strength with control.

VISIBILITY UNDER FLUOROSCOPY

All wire-guided balloons feature two radiopaque markers at the balloon shoulder, allowing for accurate placement across the stricture. The dedicated inflation port permits the introduction of contrast to further enhance visibility if desired.

TAPERED TIP

The rounded tip aids in the smooth passage of the balloon across the stricture. The wire-guided balloon is provided with a pre-loaded .035 guidewire that has a soft, flexible tip to assist in placement in tortuous or altered anatomy.

WIRE-GUIDED ESOPHAGEAL, PYLORIC, COLONIC MULTI-STAGE BALLOONS STERILE

UPN/Order Number	REF/Catalog Number	Description	Balloon Diameter (mm)	Balloon Length (cm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
MB23031	MBD-AM-BB-2	Reliant Multi-Stage Dilatation Balloon	6-7-8	5.5	230	2.8	1/Box
MB23071	MBD-BM-BB-2	Reliant Multi-Stage Dilatation Balloon	8-9-10	5.5	230	2.8	1/Box
MB23111	MBD-CM-BB-2	Reliant Multi-Stage Dilatation Balloon	10-11-12	5.5	230	2.8	1/Box
MB23151	MBD-DM-BB-2	Reliant Multi-Stage Dilatation Balloon	12-13.5-15	5.5	230	2.8	1/Box
MB23191	MBD-EM-BB-2	Reliant Multi-Stage Dilatation Balloon	15-16.5-18	5.5	230	2.8	1/Box
MB23221	MBD-FM-BB-2	Reliant Multi-Stage Dilatation Balloon	18-19-20	5.5	230	2.8	1/Box

FIXED WIRE ESOPHAGEAL MULTI-STAGE BALLOONS STERILE

UPN/Order Number	REF/Catalog Number	Description	Balloon Diameter (mm)	Balloon Length (cm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
FW24015	FWB-06/80-6/18-M	Rapide Esophageal Dilatation Balloon	6-7-8	8	180	2.8	2/Box
FW24025	FWB-08/80-6/18-M	Rapide Esophageal Dilatation Balloon	8-9-10	8	180	2.8	2/Box
FW24035	FWB-10/80-6/18-M	Rapide Esophageal Dilatation Balloon	10-11-12	8	180	2.8	2/Box
FW24045	FWB-12/80-6/18-M	Rapide Esophageal Dilatation Balloon	12-13.5-15	8	180	2.8	2/Box
FW24055	FWB-15/80-6/18-M	Rapide Esophageal Dilatation Balloon	15-16.5-18	8	180	2.8	2/Box
FW24065	FWB-18/80-6/18-M	Rapide Esophageal Dilatation Balloon	18-19-20	8	180	2.8	2/Box



ABC 3-STEP
DILATION BALLOONS

For a simple, affordable solution to balloon dilation, choose the ABC™ Dilation Balloon. Manufactured from a high-quality, non-compliant nylon material, the ABC Balloon offers the radial strength required to accurately and consistently dilate strictures.

The balloon can be inflated in 3 steps to provide gradual dilation of the stricture with a single balloon. It offers an economical solution versus typical radial expansion balloons.

The balloon is available as a flexible fixed wire esophageal 8 cm balloon or a fixed wire esophageal, pyloric and colonic 5.5 cm balloon, both featuring rapid inflation/deflation performance.



KEY BENEFITS

CHOICE OF MATERIAL
The durable nylon material delivers strength with accurate sizing. The non-compliant material provides the radial force required for consistent dilation. Inflated in 3 steps, the ABC balloon provides an affordable alternative to radial expanding balloons.

ROUNDED PROXIMAL SHOULDER
The combination of a rounded proximal shoulder and high-quality transparent nylon material facilitates visualization of the stricture during dilation.

BALLOON LENGTH
Choose from a balloon length of 5.5 cm, which can be used as a universal device to dilate esophageal, pyloric or colonic strictures, or an 8 cm balloon depending on clinical need.

ABC ESOPHAGEAL, PYLORIC, COLONIC DILATION BALLOON STERILE

UPN/Order Number	REF/Catalog Number	Description	Balloon Diameter (mm)	Balloon Length (cm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
FW24165	FWB-06/55-6/23-S	ABC Dilatation Balloon	6	5.5	230	2.8	2/Box
FW24175	FWB-08/55-6/23-S	ABC Dilatation Balloon	8	5.5	230	2.8	2/Box
FW24185	FWB-10/55-6/23-S	ABC Dilatation Balloon	10	5.5	230	2.8	2/Box
FW24195	FWB-12/55-6/23-S	ABC Dilatation Balloon	12	5.5	230	2.8	2/Box
FW24205	FWB-15/55-6/23-S	ABC Dilatation Balloon	15	5.5	230	2.8	2/Box
FW24215	FWB-18/55-6/23-S	ABC Dilatation Balloon	18	5.5	230	2.8	2/Box
FW24225	FWB-20/55-6/23-S	ABC Dilatation Balloon	20	5.5	230	2.8	2/Box

ABC ESOPHAGEAL DILATION BALLOON STERILE

UPN/Order Number	REF/Catalog Number	Description	Balloon Diameter (mm)	Balloon Length (cm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
FW24075	FWB-06/80-6/18-S	ABC Esophageal Dilatation Balloon	6	8	180	2.8	2/Box
FW24085	FWB-08/80-6/18-S	ABC Esophageal Dilatation Balloon	8	8	180	2.8	2/Box
FW24095	FWB-10/80-6/18-S	ABC Esophageal Dilatation Balloon	10	8	180	2.8	2/Box
FW24105	FWB-12/80-6/18-S	ABC Esophageal Dilatation Balloon	12	8	180	2.8	2/Box
FW24125	FWB-15/80-6/18-S	ABC Esophageal Dilatation Balloon	15	8	180	2.8	2/Box
FW24145	FWB-18/80-6/18-S	ABC Esophageal Dilatation Balloon	18	8	180	2.8	2/Box
FW24155	FWB-20/80-6/18-S	ABC Esophageal Dilatation Balloon	20	8	180	2.8	2/Box



SIMPLE SOLUTION BALLOON


INFLATION DEVICE

INFLATION DEVICE STERILE

UPN/Order Number	REF/Catalog Number	Description	Package Units
DB25011	BID-15-60-01	Micro-Tech Inflation Device 60 ml with Lock	1/Box
DB25012	BID-15-60-01	Micro-Tech Inflation Device 60 ml with Lock	5/Case

FREEDOM FROM REPROCESSING

SINGLE-USE GUIDEWIRE



Repeated use and reprocessing of reusable guidewires and cleaning brushes can result in worn, frayed, bent or otherwise damaged devices.

The use of disposable guidewires and disposable cleaning brushes can reduce the risk of device damage as well as the potential cross-contamination from using reprocessed accessories with reusable polyvinyl dilator systems.

The Single-Use Marked Spring Tip Guidewire and Dilator Cleaning Brush are for use with American and Savary dilators.

MARKED SPRING TIP GUIDEWIRE STERILE

UPN/Order Number	REF/Catalog Number	Description	Length (cm)	Package Units
DG25111	DG-18-08-21	Single-Use Marked Spring Tip Guide Wire for Esophageal Dilators	210	10/Box





DILATOR CLEANING BRUSH

UPN/Order Number	REF/Catalog Number	Description	Package Units
CB62411	CB10-0030/09-S	Dilator Cleaning Brush, 3 mm for reusable American and Savary Dilators	10/Box

EUS NEEDLES

AREUS FNA NEEDLE

TRIDENT FNB NEEDLE



DEPENDABLE & PRECISE

ENDOSCOPIC ULTRASOUND NEEDLES

Whether requiring an EUS solution for diagnostic and therapeutic intervention or for histological sampling, Micro-Tech Endoscopy provides high-quality, affordable clinical options to help with consistent sampling and accurate diagnosis from EUS tissue acquisition.

Adapt Nitinol Needles


In addition to the traditional line of cobalt chromium EUS needles, Micro-Tech now offers the Nitinol constructed Adapt needles. Nitinol provides enhanced maneuverability in tortuous scope positions and maintains its shape for multiple passes. The Trident FNB and Areus FNA designs are available in both materials.

Areus™ Fine Needle Aspiration

The Areus features an elongated dual-edged bevel tip to aid in precise penetration into the target area, delivering a consistent yield for cytology.

Trident™ Fine Needle Biopsy

The Trident has a unique tri-tip core needle design for maximized yield and enables consistent histological tissue acquisition.



KEY BENEFITS



STYLET ETCHING
Engraving on stylet provides tactile feedback when the stylet is close to exiting the needle shaft. This indicator may help prevent the stylet from releasing unexpectedly, resulting in contamination or spray.



RAPID SHEATH AND NEEDLE ADJUSTORS
The push button adjustors on Micro-Tech EUS needles offer quick and secure changes in needle depth and sheath length. The ergonomic handle design with smooth slide action provides precise one-handed control.



STYLET MANAGEMENT CAP
Integrated stylet cap that holds looped stylet securely, so it is easier to manage.

NEEDLE ECHOGENICITY
It is always important to maintain tip visibility for accurate guidance within the target area. The needle surface features a laser-etched V pattern that provides good visibility under ultrasound.

COBALT-CHROMIUM (COCR)
The Trident and Areus EUS needles are made of CoCr alloy which maintains its shape and sharpness while providing durability, flexibility and reliable target accessibility.

ADAPT NITINOL NEEDLES
Nitinol is a memory metal which returns to its original shape. Our new line of Adapt needles leverages this technology to allow greater flexibility in tortuous anatomy.

AREUS ENDOSCOPIC ULTRASOUND FINE NEEDLE ASPIRATION STERILE

UPN/Order Number	REF/Catalog Number	Description	Needle Size	Sheath Diameter (mm)	Minimum Channel Size (mm)	Package Units
ES34621	EUS-25-0-N	Areus FNA Needle CoCr	25ga	1.8	2.8	Each
ES34601	EUS-22-0-N	Areus FNA Needle CoCr	22ga	1.8	2.8	Each
ES34641	EUS-22-0-F	Areus Adapt FNA Needle Nitinol	22ga	1.8	2.8	Each
ES34631	EUS-19-0-F	Areus Adapt FNA Needle Nitinol	19ga	1.8	2.8	Each

TRIDENT ENDOSCOPIC ULTRASOUND FINE NEEDLE BIOPSY STERILE

UPN/Order Number	REF/Catalog Number	Description	Needle Size	Sheath Diameter (mm)	Minimum Channel Size (mm)	Package Units
ES34611	EUS-25-1-N	Trident FNB Needle CoCr	25ga	1.8	2.8	Each
ES34591	EUS-22-1-N	Trident FNB Needle CoCr	22ga	1.8	2.8	Each
ES34671	EUS-22-1-F	Trident Adapt FNB Needle Nitinol	22ga	1.8	2.8	Each
ES34661	EUS-19-1-F	Trident Adapt FNB Needle Nitinol	19ga	1.8	2.8	Each

- Working Length: 137.5 cm to 141.5 cm, adjustable
- Needle Length: 0 to 8 cm, adjustable
- Includes Vacuum Lock 20 cc syringe and one-way stopcock

HEMOSTASIS

LOCKADO HEMOSTASIS CLIP

DURACLIP HEMOSTASIS CLIP

SURECLIP HEMOSTASIS CLIP

RITECLIP HEMOSTASIS CLIP

PINPOINT SUBMUCOSAL NEEDLE

22-25G INJECTION NEEDLES


SUREFIRE MULTIPLE BAND LIGATOR SET

ENDOSCOPIC MULTI-BAND LIGATORS

IN-SIGHT HEMORRHOID BANDER

LOCK ONTO TISSUE


LOCKADO HEMOSTASIS CLIP



Hemostasis clips need to provide accuracy, reliability and precise functionality. Lockado® allows unrestricted opening and closing so it can be repositioned until the desired position is achieved.

With two sets of opposing atraumatic teeth in the arms of the clip, Lockado delivers a reliable and secure grip for consistent performance. Driven by Micro-Tech Endoscopy's proven clip delivery system, the Lockado clip provides exceptional reliability and tissue management.

- Jaw style works for multiple clinical applications
- Can be used with side-viewing endoscope
- Shorter stem makes the clip less obtrusive
- Variety of sizes: 11 mm, 16 mm and MAX 22 mm arm opening widths
- Indicated for anchoring jejunal feeding tubes



KEY BENEFITS

ATRAUMATIC TEETH
Projections on the clip arms help grip tissue and lock it into place.

RIBBED CATHETER
The lubricious ribbed catheter surrounds a torque wire that provides 1:1 response and 360° rotation without inadvertent clip deployment or helicoptering.

ROTATION GRIP
Secure grip allows for 360° bidirectional controlled rotation and placement of clip on targeted tissue.



LOCKADO HEMOSTASIS CLIP STERILE

UPN/Order Number	REF/Catalog Number	Opening Width (mm)	Sheath Diameter (mm)	Working Length (cm)	Minimum Channel Size (mm)	Packing Units
LO30861	LOCK-D-26-235-C-N	11	2.6	235	2.8	10/Box
LO30865	LOCK-D-26-235-C-N	11	2.6	235	2.8	2/Box



LOCKADO PLUS HEMOSTASIS CLIP STERILE

UPN/Order Number	REF/Catalog Number	Opening Width (mm)	Sheath Diameter (mm)	Working Length (cm)	Minimum Channel Size (mm)	Packing Units
LO30891	LOCK-F-26-235-C-N	16	2.6	235	2.8	10/Box
LO30895	LOCK-F-26-235-C-N	16	2.6	235	2.8	2/Box



LOCKADO MAX HEMOSTASIS CLIP* STERILE

UPN/Order Number	REF/Catalog Number	Opening Width (mm)	Sheath Diameter (mm)	Working Length (cm)	Minimum Channel Size (mm)	Packing Units
LO35041	LOCK-G-26-235-C-N	22	2.6	235	2.8	10/Box
LO35045	LOCK-G-26-235-C-N	22	2.6	235	2.8	2/Box



DEPLOY WITH PRECISION

DURACLIP™ HEMOSTASIS CLIP

Now part of the Micro-Tech portfolio, DuraClip™ continues to deliver the trusted performance and quality clinicians expect.

The DuraClip Hemostasis Clip delivers precise, reliable tissue approximation across a range of GI applications. Designed for control and simplicity, DuraClip can be opened and closed an unlimited number of times before deployment, helping clinicians reposition with confidence.

With responsive rotation for accurate positioning, flexible access in challenging anatomy, and intuitive deployment, DuraClip is designed to support efficient, controlled clip placement.

KEY FEATURES

CONTROLLED CLIP PLACEMENT

True 1:1 rotation with unlimited open and close capability before deployment supports precise positioning.

ACCESS CHALLENGING ANATOMY


Flexible distal end enables rotation over ERCP scope elevator and accommodates insertion and rotation in retroflexed scope.

DESIGNED FOR PROCEDURAL EFFICIENCY

Small stem profile and intuitive deployment support visualization and ease of use during clip placement.

DURACLIP HEMOSTASIS CLIP

UPN/Order Number	REF/Catalog Number	Opening Width (mm)	Sheath Diameter (mm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
RC40011	DC0235	11	2.6	235	2.8	10/Box
RC40031	DC0235W	16	2.6	235	2.8	10/Box



CLIP WITH CONFIDENCE

SURECLIP

Clips need to be reliable, accurate, and flexible enough to reposition or rotate as needed.

SureClip® achieves this by design delivering outstanding repositionability and reliable rotation prior to deployment, in various scope positions. SureClip's short stem aids placement in narrow lumen and improves visibility.

The proprietary clip design provides reliable deployment, may improve retention, and offers a choice of jaw sizes.

KEY BENEFITS

REPOSITIONING
SureClip's unique design allows the jaw to open, close, and reposition prior to deployment, supporting accurate placement. Fenestrations accommodate tissue and may enhance retention.

RELIABLE ROTATION
The rotation handle allows for precise clip orientation during tissue approximation or defect closure.

SHORT STEM
A shorter stem improves visualization of the target area, particularly when multiple clips are placed in close proximity.

SURECLIP^{MINI} LOWEST PROFILE STERILE

UPN/Order Number	REF/Catalog Number	Opening Width (mm)	Sheath Diameter (mm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
RC30415	ROCC-C-26-235-C-N	8	2.6	235	2.8	2/Box
RC30411	ROCC-C-26-235-C-N	8	2.6	235	2.8	10/Box

SURECLIP SHORTEST STEM STERILE

UPN/Order Number	REF/Catalog Number	Opening Width (mm)	Sheath Diameter (mm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
RC30445	ROCC-D-26-235-C-N	11	2.6	235	2.8	2/Box
RC30441	ROCC-D-26-235-C-N	11	2.6	235	2.8	10/Box

SURECLIP^{PLUS} GREATEST CLOSURE/COMPRESSION FORCE STERILE


UPN/Order Number	REF/Catalog Number	Opening Width (mm)	Sheath Diameter (mm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
RC30385	ROCC-F-26-235-C-R	16	2.6	235	2.8	2/Box
RC30381	ROCC-F-26-235-C-R	16	2.6	235	2.8	10/Box

SURECLIP^{MAX} WIDEST JAW OPENING* STERILE

UPN/Order Number	REF/Catalog Number	Opening Width (mm)	Sheath Diameter (mm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
RC30545	ROCC-H-26-235-C-N	17	2.6	235	2.8	2/Box
RC30541	ROCC-H-26-235-C-N	17	2.6	235	2.8	10/Box

*Jaws may not open to maximum reach until placed in contact with tissue

CHOOSE PERFORMANCE & ECONOMY
RITECLIP



Rest assured, with the RiteClip™, there is no need to repeat clip placements. Accurate initial clip positioning may reduce procedure time and the number of clips needed to achieve consistent results.

RiteClip delivers reliable rotation even in extended-scope positions. RiteClip can be repositioned as needed before deployment and, with its short stem design, is less obtrusive in the endoscopic field of view.

Featuring a comfortable 3-ring handle and flexible distal shaft, this proven delivery system blends both performance and cost savings.

KEY BENEFITS

REPOSITIONING
The RiteClip features the ability to open and close the clip without limitations. The ability to close the clip, assess, and if needed, adjust the clip placement contributes to improved outcomes.

RELIABLE ROTATION
RiteClip includes a flexible distal section that allows for accurate rotation, even in advanced scope positions. Rotation is simply achieved by rotating the 3-Ring handle.

SHORT STEM
A shorter stem makes the clip less obtrusive, improving visualization of the target area, particularly when multiple clips are placed in close proximity.


RITECLIP STERILE

UPN/Order Number	REF/Catalog Number	Nominal Jaw Opening (mm)	Minimum Working Channel (mm)	Working Length (cm)	Package Units
RT30665	ROCC-D-26-235-C	11	2.8	235	2/Box
RT30661	ROCC-D-26-235-C	11	2.8	235	10/Box

RITECLIP PLUS STERILE

UPN/Order Number	REF/Catalog Number	Nominal Jaw Opening (mm)	Minimum Working Channel (mm)	Working Length (cm)	Package Units
RT30785	ROCC-F-26-235-C	16	2.8	235	2/Box
RT30781	ROCC-F-26-235-C	16	2.8	235	10/Box

*Jaws may not open to maximum reach until placed in contact with tissue



SHARP INNOVATION

INJECTION NEEDLES

A simple, elegant solution, the 5 mm and 3 mm long lancet Injection Needles are offered in 22 and 25 gauge. Designed for single-handed operation, the innovative design limits accidental needle deployment during insertion and removal. The patented ratchet mechanism provides audible confirmation of deployment and minimizes needle bounce. A quick release button returns the needle fully into the sheath. The recessed metal end cap and firm catheter aid in tissue puncture and ensure needle performance through various scope positions.

KEY BENEFITS

PATENTED RATCHET DESIGN

Provides audible confirmation and secure locking of the deployed needle during use. As the needle is advanced, the ratchet gives tactile feedback and helps to reduce needle bounce.

TARGET TISSUE DEPTH

With the shorter needle tip, PinPoint helps to consistently target the submucosal space. The 3 mm projection is sufficient to reach the submucosa in even the thickest intestinal tissue.

PINPOINT INJECTION NEEDLE


Needle projection shortened to 3 mm which helps target the submucosal space.

SHORT BEVEL TIP

Both the 5 mm and 3 mm needles feature a shorter three-bevel tip, designed to improve accuracy. The 3 mm PinPoint needle was designed to help target the submucosal space while the 5 mm allows for greater flexibility in therapeutic use.

ONE-HANDED OPERATION

Ergonomically designed one-handed operation enables confident needle deployment and rapid needle retraction.



*Tissue depth representative of thickest tissue recorded in study, Echoendoscopic characterization of the human colon, Fernando M. Castro-Poças, Mário Dinis-Ribeiro, Tarcisio P. Araújo and Isabel Pedroto

PINPOINT™ SUBMUCOSAL STERILE

UPN/Order Number	REF/Catalog Number	Description	Needles Gauge/Length (mm)	Sheath Diameter (mm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
IN33301	IN12-223232302	PinPoint Injection Needle Lower	22/3	2.3	230	2.8	10/Box
IN33331	IN12-223232302	PinPoint Injection Needle Lower	25/3	2.3	230	2.8	10/Box

INJECTION NEEDLES STERILE

UPN/Order Number	REF/Catalog Number	Description	Needles Gauge/Length (mm)	Sheath Diameter (mm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
IN33131	IN12-225231802	Injection Needle Upper	22/5	2.3	180	2.8	10/Box
IN33221	IN12-225231802	Injection Needle Upper	25/5	2.3	180	2.8	10/Box
IN33151	IN12-225232302	Injection Needle Lower	22/5	2.3	230	2.8	10/Box
IN33241	IN12-225232302	Injection Needle Lower	25/5	2.3	230	2.8	10/Box



SIMPLE SETUP. CONFIDENT DEPLOYMENT.

SUREFIRE™ MULTIPLE BAND LIGATOR SET

Band ligation setup shouldn't be complicated. It should be intuitive, dependable and ready to use without extra steps.

The SureFire™ Multiple Band Ligator Set is designed to simplify setup with its preloaded, over-the-scope design that reduces setup steps while keeping the working channel free. Audible and tactile feedback, paired with clear band numbering, provide added confidence during deployment and support a reliable, streamlined banding process.



KEY BENEFITS

PRELOADED, OVER-THE-SCOPE DESIGN

Simplifies setup while keeping the working channel free to support a more efficient workflow.

CLEAR VISUAL FIELD

Designed to maintain a clear procedural view during band deployment, without obstructing the primary viewing area.

AUDIBLE & TACTILE FEEDBACK

Provides audible and tactile confirmation, allowing you to hear and feel each band release.

BAND NUMBERING INDICATION

Clear numbering on the handle helps track each band in use, supporting confident, sequential deployment.

SUREFIRE™ MULTIPLE BAND LIGATOR SET

Order Number	REF/Catalog Number	Description	Working Length (cm)	Package Units
ML32931	MBLS-6-NL-S	Endoscopic Multi-Band Ligator Set, 6 Synthetic Rubber Bands, compatible with 9.5-11mm scope diameter	110	1/Box

ENDOSCOPIC MULTIPLE BAND LIGATORS



Endoscopic banding has been demonstrated as an effective and simple treatment for internal hemorrhoids and varices. Hemorrhoid banding is readily performed in an ambulatory setting and can result in less pain and have a shorter recovery period than a surgical hemorrhoidectomy.

The Multi-Band Ligator Sets from Micro-Tech offer cost-effective solutions for the treatment of hemorrhoids and varices during routine endoscopy procedures.

KEY BENEFITS

EXCELLENT VISIBILITY

A clear ligating barrel with rear-mounted bands optimizes the endoscopic view while maintaining suction.

BAND DESIGN

The square band design provides consistent shape and reliable retention. The hemorrhoid set comes with four pre-loaded latex-free bands for routine colonoscopies. The upper ligator comes pre-loaded with six latex-free bands and can be used in upper and lower endoscopic procedures.

USER-FRIENDLY SYSTEM

The handle unit comes with the pull wire pre-loaded to simplify setup. The alert band notifies the operator when only one band remains.

MULTI-BAND LIGATOR STERILE ②

UPN/Order Number	REF/Catalog Number	Description	Minimum Working Channel (mm)	Working Length (cm)	Package Units
HL32161	HMBLS-4F-NL	Endoscopic Hemorrhoid Multi-Band Ligator Set, 4 Latex-Free Bands, compatible with 9.4-13 mm scope diameter	2.8	220	1/Box
ML32651	MBLS-6F-NL	Endoscopic Multi-Band Ligator Set, 6 Latex-Free Bands, compatible with 9.4-13 mm scope diameter	2.8	145	1/Box

IN-OFFICE LIGATION

IN-SIGHT BANDER



Hemorrhoid banding has been demonstrated as an effective and simple treatment for internal hemorrhoids, and can be readily performed in an office or ambulatory setting. Banding can also result in less pain with a shorter recovery period than a surgical hemorrhoidectomy.

The In-Sight™ Multi-Band Ligator is the only hand-held ligation device with its own light source – providing illumination, suction, direct visualization and ligation capabilities in a single-use device.

- Preloaded with 4 ligator bands
- Comfortable and ergonomic design for easy ligation
- Unique LED system maximizes visualization safely and efficiently



IN-SIGHT MULTI-BAND LIGATOR

UPN/Order Number	REF/Catalog Number	Description	Package Units
QL30011	MDHL100000A	In-Sight Multi-Band Ligator	10/Box

INFECTION PREVENTION

ENDOSCOPE VALVES AND ACCESSORIES

BIOPSY VALVES

VALVES AND ACCESSORIES

DISPOSABLE TUBING

BITE BLOCKS

DISTAL TIP PROTECTORS

CLEANING BRUSHES

GI PROCEDURE KITS

ENHANCED INFECTION PREVENTION

SINGLE-USE ENDOSCOPE VALVES

Reusable medical devices, like endoscope valves, can pose a risk of infection if not adequately reprocessed between patient uses. The manual cleaning and reprocessing required to ensure complete sterilization can be labor-intensive and time-consuming.

Our disposable valves offer an efficient alternative to reusable valves, eliminating the need for extensive cleaning and reprocessing. Using our disposable valves can help improve procedural efficiency and reduce the risk of cross-contamination, enhancing your infection prevention practices.



NON-LUBRICATED SINGLE-USE ENDOSCOPE VALVES STERILE

UPN/Order Number	REF/Catalog Number	Description	Package Units
QL71151	AW01-OLY3	Air/Water & Suction, Biopsy Valve Kit, Olympus®, 3-Piece Kit	50/Box
QL71101	AW01-OLY4	Air/Water & Suction, Biopsy Valve, Aux Water Kit, Olympus®, 4-Piece Kit,	50/Box
QL71112	AW01-OLY5	Air/Water & Suction, Biopsy Valve, Aux Water, Cleaning Adapter Kit, Olympus®, 5-Piece Kit, (4 Boxes of 10)	40/Case

JACKET DISPOSABLE BIOPSY VALVE



The Jacket™ Disposable Biopsy Valve is used to cover the opening to the biopsy/suction channel of gastrointestinal endoscopes.

Its unique design offers improved valve performance by minimizing leakage and helping to preserve insufflation to maintain adequate visualization.

- Double valve system to minimize leakage and prevent spraying
- Improved seal with three points of contact
- Internal funnel seal hugs and flexes with inserted accessories without gaping
- Compatible with Olympus®, Pentax®, FujiFilm™ GI endoscopes



JACKET DISPOSABLE BIOPSY VALVE

UPN/Order Number	REF/Catalog Number	Description	Package Units
QL70021	MD-A-BV02	Jacket Disposable Biopsy Valve, Olympus®/Fujifilm™, Non-Sterile (Blue)	100/Box
QL70051	MD-B-BV02	Jacket Disposable Biopsy Valve, Pentax® Medical, Non-Sterile (Clear)	100/Box

STANDARD DISPOSABLE BIOPSY VALVES

STERILE

UPN/Order Number	REF/Catalog Number	Description	Package Units
QL72191	BW01-OLYB	Disposable Biopsy Valve, Olympus®/Fujifilm™, Sterile (Light Blue)	100/Box
QL72201	BW01-PENB	Disposable Biopsy Valve, Pentax® Medical, Sterile (Black)	100/Box

ENHANCED SAFETY AND EFFICIENCY

SINGLE-USE ENDOSCOPE VALVES AND ACCESSORIES



Reusable medical devices, like endoscope valves, can present numerous challenges for patients and clinicians, as they can become contaminated, posing an infection risk if not adequately reprocessed between patient uses. Additionally, the manual cleaning and reprocessing of these devices can be time-consuming, involving multiple steps to ensure complete sterilization and traceability to a specific endoscope.

Designed to address these challenges, our disposable valves and accessories are a reliable alternative to reusable devices that can help streamline procedures and enhance infection prevention practices by eliminating the need for timely cleaning, reprocessing, and traceability procedures and reduce the potential risk of infection from cross-contamination.

COMPREHENSIVE OFFERINGS TO MEET YOUR NEEDS

Our single-use valves and accessories are compatible with Olympus®, Fujifilm®, and Pentax® Endoscopes, ensuring seamless integration with your existing equipment. Available in 3-, 4-, or 5-piece kits and single-pack configurations, our comprehensive offerings feature air/water valves, suction and biopsy valves, auxiliary water connectors, and cleaning adapters, providing the right solution to meet your needs.

OLYMPUS® ENDOSCOPES



FUJIFILM® 700/800 SERIES ENDOSCOPES



PENTAX® ENDOSCOPES



LUBRICATED SINGLE-USE ENDOSCOPE VALVES & ACCESSORIES

UPN/Order Number	REF/Catalog Number	Description	Package Units
OLYMPUS® ENDOSCOPES			
QL72111	BW01-OLY5	5-Piece Kit (Air/Water, Suction, & Biopsy Valves w/ Auxiliary Water Connector & Cleaning Adapter)	50/Box
QL72101	BW01-OLY4	4-Piece Kit (Air/Water, Suction, & Biopsy Valves w/ Auxiliary Water Connector)	50/Box
QL72151	BW01-OLY3	3-Piece Kit (Air/Water, Suction, & Biopsy Valves)	50/Box
QL72191	BW01-OLYB	Single-Pack Biopsy Valve (Olympus®/Fujifilm® Compatible)	100/Box
QL72171	BW01-OLYA	Single-Pack Auxiliary Water Connector	100/Box
QL72231	BW01-OLYC	Single-Pack Cleaning Adapter*	100/Box
FUJIFILM® 700/800 ENDOSCOPES			
QL72131	BW01-FUJ5-F7	5-Piece Kit (Air/Water, Suction, & Biopsy Valves w/ Auxiliary Water Connector & Cleaning Adapter)	50/Box
QL72281	BW01-FUJ4-F7	4-Piece Kit (Air/Water, Suction, & Biopsy Valves w/ Auxiliary Water Connector)	50/Box
QL72301	BW01-FUJ3-F7	3-Piece Kit (Air/Water, Suction, & Biopsy Valves)	50/Box
QL72191	BW01-FUJB	Single-Pack Biopsy Valve (Olympus®/Fujifilm® Compatible)	100/Box
QL72211	BW01-FUJA	Single-Pack Auxiliary Water Connector	100/Box
QL72261	BW01-FUJC-F7	Single-Pack Cleaning Adapter*	100/Box
PENTAX® ENDOSCOPES			
QL72121	BW01-PEN4	4-Piece Kit (Air/Water, Suction, & Biopsy Valves w/ Auxiliary Water Connector)	50/Box
QL72161	BW01-PEN4	3-Piece Kit (Air/Water, Suction, & Biopsy Valves)	50/Box
QL72201	BW01-PENB	Single-Pack Biopsy Valve	100/Box
QL72181	BW01-PENA	Single-Pack Auxiliary Water Connector	100/Box

*Single-pack cleaning adapter is non-sterile.

All trademarks are the property of their respective owners.

ENHANCED EFFICIENCY AND PROTECTION

DISPOSABLE GI TUBING SOLUTIONS



IRRIGATION TUBING

Infection prevention is crucial for ensuring patient safety during endoscopic procedures. Our Disposable Tubing Solutions are specifically designed to minimize the risk of cross-contamination while enhancing procedural efficiency.

With reliable fluid delivery and optimal visibility, our tubing solutions offer healthcare professionals exceptional and easy-to-use tools that can help enhance infection prevention practices.



IRRIGATION TUBING

- 24-hour multi-patient use without need for reprocessing or tracking.
- Supports improved patient outcomes by providing a clear field of view in procedures.
- Features a Luer-lock valve designed to reduce contamination risks while maintaining optimal fluid pressure.



AIR/WATER TUBING

- 24-hour multi-patient use without need for reprocessing or tracking.
- Designed to facilitate consistent air and water flow, enhancing the clarity and efficiency of endoscopic procedures.
- Supports room air or CO2 for insufflation.



AUXILIARY WATER CONNECTORS

- Provides secure, leak-proof connections, allowing for consistent flow of fluids, enhancing the reliability of procedures.
- Easy attachment/detachment, facilitating rapid setup and teardown, minimizing downtime during procedures.

KEY BENEFITS

DISPOSABLE DESIGN
Provides 24-hour multi-patient use, minimizing the risk of cross-contamination by eliminating the need for reprocessing.

EASE OF USE
Streamlined design that allows for quick setup and easy disposal, enhancing procedural efficiency.

VERSATILE COMPATIBILITY
Connects seamlessly with sterile water bottles*, irrigation pumps and supports room air and CO2 for insufflation.

FLEXIBLE MATERIALS
Designed to withstand bends without compromising integrity, providing durability and reliability during use.

OPTIMAL FLOW
Provides consistent and optimal fluid flow giving way to clear visibility during procedures.

Order Number	Reference Number	Description	Package Units
AIR/WATER TUBING			
QL76101	BW02-OLY	Air/Water Bottle Cap Tubing Olympus® Endoscopes	10/Box
QL76111	BW02-OLYC02	Air/Water Bottle Cap CO2 Tubing Olympus® Endoscopes	15/Box
IRRIGATION TUBING			
QL76121	BW02-IRRTUB	Aux Water Channel Tubing Set Daily (compatible with Olympus® OFP-2 Irrigation Pumps)	15/Box
QL76151	BW02-IRRTUBSL	Aux Water Channel Tubing Set Daily (compatible with Olympus® AFU-100, ERBE EIP™ 2 , Endogator® EGP-100, Endo Stratus® Irrigation Pumps)	15/Box
AUXILIARY WATER CONNECTORS			
QL72211	BW01-FUJA	Single-Pack Auxiliary Water Connector, Fujifilm® Endoscopes	100/Box
QL72171	BW01-OLYA	Single-Pack Auxiliary Water Connector, Olympus® Endoscopes	100/Box
QL72181	BW01-PENA	Single-Pack Auxiliary Water Connector, Pentax® Endoscopes	100/Box
CO2 SOURCE TUBING WITH INLINE FILTER			
QL76311	BW02-CO2TUBOLY	CO2 Source Tubing, Olympus®	15/Box
QL76321	BW02-CO2TUBAOLY	CO2 Source Tubing Adapter, Olympus®	15/Box
CO2 SOURCE TUBING ADAPTER, REUSABLE			
QL76211	BW02-CO2CONOLY	CO2 Source Tubing Adapter Olympus® UCR	1/Box
QL76221	BW02-CO2CONFJ	CO2 Source Tubing Adapter Fuji® GW100	1/Box

All trademarks are the property of their respective owners.
*Compatible with B. Braun® and Baxter® sterile water bottles.

PREMIUM MATERIAL AND ENHANCED DESIGN

BREATHE-EZ BITE BLOCK



Micro-Tech Bite Blocks offer comfort and protection during upper endoscopic procedures. The new Breathe-EZ™ Bite Block was designed with a special, softer bite plate that helps to keep the lower jaw forward and the airway open.

- Unique design that may assist in improving oxygenation during endoscopic procedures
- Available with adjustable elastic strap
- Softer overlay that helps protect patient's teeth

BREATHE-EZ BITE BLOCK

UPN/Order Number	REF/Catalog Number	Description	Package Units
QL26031	MDBB1AY003A	Breathe-EZ Bite Block, Adult, 60Fr, with elastic strap	50/Box

STANDARD BITE BLOCKS

Micro-Tech Bite Blocks offer comfort and protection during upper endoscopic procedures. Available in 38Fr, 60Fr and 60+Fr sizes, with or without dental rim.

UPN/Order Number	REF/Catalog Number	Description	Package Units
BB60011	AC01-103.A	Bite Block, Adult, 60Fr, 26 mm x 21 mm opening, with strap	100/Box
BB60021	AC01-103.P	Bite Block, Pediatric, 38Fr, 20mm x 13mm opening, with strap	100/Box
BB60041	AC01-108.M	Bite Block, Adult w/o dental rim. 60+Fr, 26 mm x 23 mm opening, with strap	100/Box
BB60051	AC01-109.M	Bite Block, Adult, 60+Fr, 26 mm x 23 mm opening, with strap	100/Box

MINIMIZE SCOPE DAMAGE

ENDOSCOPE DISTAL TIP PROTECTORS



Protecting your valuable equipment is critical to maintaining high-quality patient care and operational efficiency. However, the risk of damage to endoscope tips during handling, storage, or transportation can lead to costly repairs and cause unnecessary downtime.

Our Distal Tip Protectors offer a cost-effective, easy-to-use solution to these problems. Designed for single use, our scope protectors safeguard the delicate distal end of your scopes, so you can continue to deliver exceptional patient care.

Available in two sizes:

- Large Endoscope Distal Tip Protector (blue) is designed for endoscopes with 7 mm to 15 mm in diameter.
- Small Endoscope Distal Tip Protector (red) is designed for endoscopes with 2 mm to 10 mm in diameter.

KEY FEATURES

COMPREHENSIVE PROTECTION

Designed to protect the distal end of endoscopes during handling, storage and transport, minimizing the risk of costly accidental damage.

USER FRIENDLY

Quick and simple application enables seamless integration into your current workflow.

ENDOSCOPE DISTAL TIP PROTECTORS

UPN/Order Number	REF/Catalog Number	Description	Compatible scope OD (mm)	Package Units
QL78111	BW03-DTPL	Distal Tip Protector Large (blue)	7 mm to 15 mm	100/Box
QL78101	BW03-DTPS	Distal Tip Protector Small (red)	2 mm to 10 mm	100/Box

PROTECT, COLLECT, CLEAN

INFECTION PROTECTION



CATCHALL SCOPE TRANSPORT PAD

The CatchAll™ Transport Pad offers compliance and protection while transporting contaminated scopes.

UPN/Order Number	REF/Catalog Number	Description	Package Units
CP62131	CTP-02025	CatchAll Scope Transport Pad, Blue	25/Box







CLEANING BRUSHES

Micro-Tech's Cleaning Brushes combine performance and compliance. Available in combination kits, the brush range features 6/10 mm nylon brush heads designed for effective, gentle cleaning.

UPN/Order Number	REF/Catalog Number	Description	Package Units
QL22011	MDBR1A0504A	Single-ended 6 mm channel brush	50/Box
QL22021	MDBR1B0016A	Double-ended 6 mm channel brush	50/Box
QL22031	MDBR1B0015A	Single-ended 6 mm channel and attached 10 mm valve brush	50/Box
QL22041	MDBR1B0005A	Double-ended 5 mm-10 mm valve brush	50/Box
QL22051	MDBR1B0013A	Single-ended 6 mm channel brush + 5 mm-10 mm valve brush	50/Box
QL22061	MDBR1B0014A	Double-ended 6 mm channel brush + 5 mm-10 mm valve brush	50/Box

CONVENIENT COMPLIANCE
GI PROCEDURE KITS



Procedure kits provide the essential supplies needed for most cases in one convenient package. Using preconfigured accessory kits helps to standardize processes, reduce product handling, and simplify supply ordering.

Choose four popular kits that provide the items required for nearly every endoscopic procedure. By bundling routine items into standardized kits, the number of touchpoints, and the potential for cross-contamination may be reduced.

	Catch-All Transport Pad	Biohazard Label	Bowl	Sponge (w/ detergent)	Gauze (10)	Lubricating Gel (1.4 oz)	Cleaning Brush*	Suction Tubing*
Compliance + Kit	•	•	•		•	•		
Basic Kit	•	•	•	•	•	•		
Basic + Kit	•	•	•	•	•	•	•	
Set-Up Kit	•	•	•	•	•	•	•	•

*Cleaning Brush (Single-ended 6 mm channel and attached 10 mm valve brush); Suction Tubing (3/16" x 10')

UPN/Order Number	REF/Catalog Number	Description	Package Units
QL22561	KT1A0006A	Compliance + Endoscopy Kit	30/Box
QL22511	MDKT1A0001A	Basic Endoscopy Kit	12/Box
QL22551	MDKT1A0005A	Basic + Endoscopy Kit	12/Box
QL22531	MDKT1A0003A	Set-Up Endoscopy Kit	12/Box
QL22711	MDJL1A0001A	Lubricating Gel, 1.4 oz Pouch	100/Box
QL22811	MDSB1B0001A	Endoscope Enzymatic Cleaning Sponge	50/Box

NITINOL STENTS

SEGMENTED ESOPHAGEAL STENTS

FULLY & PARTIALLY COVERED ESOPHAGEAL STENTS

PRECISE PLACEMENT


SEGMENTED ESOPHAGEAL STENTS

Precise stent placement and maintaining position throughout the deployment process are essential. Micro-Tech Endoscopy has over two decades of experience making high-quality, affordable esophageal stent options. Our segmented esophageal stent conforms to bypass esophageal stenoses and allows greater flexibility, which may be more comfortable for patients.

The deployment system has a pliable and tapered tip to facilitate placement through tight strictures while the pre-loaded delivery catheter and sliding handle mechanism provide smooth delivery of the stent to the targeted site. Our line of esophageal stents combine flexibility and consistent radial force to help with stricture resolution.



KEY BENEFITS



RADIOPAQUE MARKER
Although the entire stent is visible under X-ray fluoroscopy, each stent has several radiopaque markers to assist with precise stent placement.

MIGRATION RESISTANCE
The flared ends along with the unique segmented construction may reduce migration within the esophageal lumen.

INDEPENDENT MOVEMENT
The segmented modular stent allows independent movement between each section which maximizes conformance to curved or tortuous anatomy.

REPOSITIONING SUTURES
Sutures are located on the distal and proximal ends of the stent to allow repositioning or removal immediately post stent placement.



SEGMENTED ESOPHAGEAL STENT STERILE ②

UPN/ Order Number	REF/Catalog Number	Description	Diameter of Stent (mm)	Diameter of Flares (mm)	Stent Length (cm)	Catheter Outer Diameter (mm)	Working Length (cm)
ET80510	NST61-224-18.060	Segmented Esophageal Stent*	18	24	6	6.7	65
ET80520	NST61-224-18.080	Segmented Esophageal Stent*	18	24	8	6.7	65
ET80530	NST61-224-18.100	Segmented Esophageal Stent*	18	24	10	6.7	65
ET80540	NST61-224-18.120	Segmented Esophageal Stent*	18	24	12	6.7	65
ET80550	NST61-224-18.140	Segmented Esophageal Stent*	18	24	14	6.7	65
ET80560	NST61-224-20.060	Segmented Esophageal Stent*	20	26	6	6.7	65
ET80570	NST61-224-20.080	Segmented Esophageal Stent*	20	26	8	6.7	65
ET80580	NST61-224-20.100	Segmented Esophageal Stent*	20	26	10	6.7	65
ET80590	NST61-224-20.120	Segmented Esophageal Stent*	20	26	12	6.7	65
ET80600	NST61-224-20.140	Segmented Esophageal Stent*	20	26	14	6.7	65
ET80610	NST61-224-22.060	Segmented Esophageal Stent*	22	28	6	6.7	65
ET80620	NST61-224-22.080	Segmented Esophageal Stent*	22	28	8	6.7	65
ET80630	NST61-224-22.100	Segmented Esophageal Stent*	22	28	10	6.7	65
ET80640	NST61-224-22.120	Segmented Esophageal Stent*	22	28	12	6.7	65
ET80650	NST61-224-22.140	Segmented Esophageal Stent*	22	28	14	6.7	65


Note: The inner tube contains a central lumen that accommodates a .038 in (.97 mm) guide wire. *Non-Stock Item

PRECISE PLACEMENT

FULLY COVERED & PARTIALLY COVERED ESOPHAGEAL STENTS

Precise stent placement and maintaining position throughout the deployment process are essential. Micro-Tech Endoscopy has over two decades of experience making high-quality, affordable esophageal stent options.

The deployment system has a pliable and tapered tip to facilitate placement through tight strictures, while the pre-loaded delivery catheter and sliding handle mechanism provide smooth delivery of the stent to the targeted site.



KEY BENEFITS

RADIOPAQUE MARKERS

Although the entire stent is visible under X-ray fluoroscopy, each stent has several radiopaque markers to assist with precise stent placement.

MANUFACTURED FROM ONE WIRE

Our fully and partially covered esophageal stents are braided from a single Nitinol wire, reducing the risk of fracture at the connection points.

MIGRATION RESISTANCE

The flared ends may reduce migration within the esophageal lumen.

PARTIALLY COVERED ESOPHAGEAL STENT STERILE



UPN/Order Number	REF/Catalog Number	Description	Diameter of Stent (mm)	Diameter of Flares (mm)	Stent Length (cm)	Catheter Outer Diameter (mm)	Working Length (cm)
ET80010	NST01-222-18.060	Partially Covered Esophageal Stent*	18	23	6	6.7	65
ET80020	NST01-222-18.080	Partially Covered Esophageal Stent*	18	23	8	6.7	65
ET80030	NST01-222-18.100	Partially Covered Esophageal Stent*	18	23	10	6.7	65
ET80040	NST01-222-18.120	Partially Covered Esophageal Stent*	18	23	12	6.7	65
ET80050	NST01-222-18.140	Partially Covered Esophageal Stent*	18	23	14	6.7	65
ET80060	NST01-222-18.150	Partially Covered Esophageal Stent*	18	23	15	6.7	65
ET80070	NST01-222-20.060	Partially Covered Esophageal Stent*	20	25	6	6.7	65
ET80080	NST01-222-20.080	Partially Covered Esophageal Stent*	20	25	8	6.7	65
ET80090	NST01-222-20.100	Partially Covered Esophageal Stent*	20	25	10	6.7	65
ET80100	NST01-222-20.120	Partially Covered Esophageal Stent*	20	25	12	6.7	65
ET80110	NST01-222-20.140	Partially Covered Esophageal Stent*	20	25	14	6.7	65
ET80120	NST01-222-20.150	Partially Covered Esophageal Stent*	20	25	15	6.7	65
ET80130	NST01-222-22.060	Partially Covered Esophageal Stent*	22	27	6	6.7	65
ET80140	NST01-222-22.080	Partially Covered Esophageal Stent*	22	27	8	6.7	65
ET80150	NST01-222-22.100	Partially Covered Esophageal Stent*	22	27	10	6.7	65
ET80160	NST01-222-22.120	Partially Covered Esophageal Stent*	22	27	12	6.7	65
ET80170	NST01-222-22.140	Partially Covered Esophageal Stent*	22	27	14	6.7	65
ET80180	NST01-222-22.150	Partially Covered Esophageal Stent*	22	27	15	6.7	65

FULLY COVERED ESOPHAGEAL STENT STERILE



UPN/Order Number	REF/Catalog Number	Description	Diameter of Stent (mm)	Diameter of Flares (mm)	Stent Length (cm)	Catheter Outer Diameter (mm)	Working Length (cm)
ET80190	NST01-224-18.060	Fully Covered Esophageal Stent*	18	23	6	6.7	65
ET80200	NST01-224-18.080	Fully Covered Esophageal Stent*	18	23	8	6.7	65
ET80210	NST01-224-18.100	Fully Covered Esophageal Stent*	18	23	10	6.7	65
ET80220	NST01-224-18.120	Fully Covered Esophageal Stent*	18	23	12	6.7	65
ET80230	NST01-224-18.140	Fully Covered Esophageal Stent*	18	23	14	6.7	65
ET80240	NST01-224-18.150	Fully Covered Esophageal Stent*	18	23	15	6.7	65
ET80250	NST01-224-20.060	Fully Covered Esophageal Stent*	20	25	6	6.7	65
ET80260	NST01-224-20.080	Fully Covered Esophageal Stent*	20	25	8	6.7	65
ET80270	NST01-224-20.100	Fully Covered Esophageal Stent*	20	25	10	6.7	65
ET80280	NST01-224-20.120	Fully Covered Esophageal Stent*	20	25	12	6.7	65
ET80290	NST01-224-20.140	Fully Covered Esophageal Stent*	20	25	14	6.7	65
ET80300	NST01-224-20.150	Fully Covered Esophageal Stent*	20	25	15	6.7	65
ET80310	NST01-224-22.060	Fully Covered Esophageal Stent*	22	27	6	6.7	65
ET80320	NST01-224-22.080	Fully Covered Esophageal Stent*	22	27	8	6.7	65
ET80330	NST01-224-22.100	Fully Covered Esophageal Stent*	22	27	10	6.7	65
ET80340	NST01-224-22.120	Fully Covered Esophageal Stent*	22	27	12	6.7	65
ET80350	NST01-224-22.140	Fully Covered Esophageal Stent*	22	27	14	6.7	65

Note: The inner tube contains a central lumen that accommodates a .038 in (.97 mm) guide wire. *Non-Stock Item

EMR-ESD

LESIONHUNTER

EDGEHOG

ADVANTIS

SURETRAC

BLUEBOOST

GOLDKNIFE

ENSURE

DUAL ACTION CLOSURE

LOCKADO MAX

SPECIMEN TRAY


DEVICE ORGANIZER

PRECISION CUT

LESIONHUNTER COLD SNARE

LesionHunter™ Cold Snare features a rotatable, thin Nitinol Wire which maintains an oval shape cut after cut. With this cold snare, you get a clean and precise cut.

- Dedicated Cold Snare – Specifically designed for resecting polyps using a cold snare technique
- Thin Nitinol Wire – This wire material provides durability in the shape and allows for a 0.18 mm thin wire for an even cleaner cut
- Coiled Sheath – The coiled metal catheter creates greater stability for the guillotine-like cutting action of the wire
- Comfortable Ergonomic Handles – Like other Micro-Tech snares, this snare has a three-ring comfortable grip with shortened throw for an effortless cut
- Reduced Fly-Away Tissue – The thin wire may reduce the risk of generating a "fly-away" specimen, thereby making it easier to find the tissue at the end of the resection



LESIONHUNTER STERILE

UPN/Order Number	REF/Catalog Number	Description	Length (cm)	Outer Diameter (mm)	Package Units
CS50051	CS2-21023231	LesionHunter Cold Snare, Oval 10 mm wide, Nitinol Wire	230	2.3	10/Box
CS50061	CS2-21523231	LesionHunter Cold Snare, Oval 15 mm wide, Nitinol Wire	230	2.3	10/Box
CS50071	CS2-22023231	LesionHunter Cold Snare, Oval 20 mm wide, Nitinol Wire	230	2.3	10/Box

PRECISION RESECTION

EDGEHOG HOT SNARE

EdgeHog™ is a rotatable hot snare featuring a flat monofilament wire and ultra-stiff construction. The device is designed to provide controlled, multi-stage deployment, allowing the snare to be positioned at different loop diameters during polypectomy and EMR procedures. Its rounded silhouette and wire configuration are intended to support positioning, manipulation, and tissue capture during endoscopic resection.

TISSUE ENGAGEMENT

The flat monofilament wire is designed to engage tissue during snare closure. The wire is oriented perpendicular to the cutting plane, while the proximal V-shaped wire configuration is intended to help maintain loop shape and stability throughout deployment and closure. These design features are intended to support tissue capture and snare positioning during resection procedures.

FLAT MONOFILAMENT WIRE

The flat wire is designed to provide consistent energy delivery during electrosurgical resection. The wire profile and configuration are intended to support controlled tissue resection and visualization of the treatment area during endoscopic procedures.



EDGEHOG HOT SNARE

UPN/Order Number	REF/Catalog Number	Description	Loop Diameter	Working Length	Package Units
PS51181	PS51181	EdgeHog Hot Snare	20 mm	230 cm	10/Box


PRECISION RESECTION

ADVANTIS FTS



The ADVANTIS snare is the newest addition to the Micro-Tech snare portfolio. This snare features a firm braided wire and a domed tip providing:

- Feature-driven design for polypectomy and EMR
- Stiff 1x7 braided wire structure offers strength to meaningfully manipulate tissue in loop closure
- Domed snare tip for precise, procedure-specific techniques
- 15 mm or 20 mm loop diameter options yield en bloc or piecemeal resection capabilities
- 0.42 mm wire diameter provides strength to snare loop structure and integrity



ADVANTIS FTS

UPN/Order Number	REF/Catalog Number	Description	Loop diameter (mm)	Working Length (cm)	Package Units
PS52051	NRPFS05-01523230	Advantis Small FTS Oval	15	230	10/Box
PS52061	NRPFS05-02023230	Advantis Medium FTS Oval	20	230	10/Box



VISUALIZED DISECTION

SURETRAC TRACTION SYSTEM

SureTrac™ retracts tissue to assist in tissue exposure, resection, and removal within the stomach and colon. The SureTrac system features a preloaded elastic silicone band on a primary clip and a secondary clip to anchor on tissue. The superelastic silicone band is designed to retract and release elastic potential energy for dynamic traction; this provides enhanced visualization and enables tissue management while dissecting submucosal lesions.

Each kit includes:

- A Primary Device, SureClip® preloaded with a silicone band
- A Secondary Device with the same size clip as the Primary Device but available in either SureClip® or Lockado® clip design.

SURETRAC KIT STERILE

UPN/Order Number	REF/Catalog Number	Description	Clip Size	Clipping Device Style	Minimum Channel Size	Package Units
TS63041	ETS00004	SureTrac Elastic Traction Device	11 mm	SureClip	2.8 mm	1/Box
TS63061	ETS00006	SureTrac Elastic Traction Device	16 mm	SureClip	2.8 mm	1/Box
TS63101	ETS00004	SureTrac Elastic Traction Device	11 mm	Lockado	2.8 mm	1/Box
TS63121	ETS00006	SureTrac Elastic Traction Device	16 mm	Lockado	2.8 mm	1/Box

CLEAR DIFFERENTIATION

BLUEBOOST SUBMUCOSAL INJECTABLE SOLUTION



BlueBoost™ is a lifting agent that lasts. Once injected into the submucosal space, BlueBoost adds efficiency to GI procedures by separating the layers of the intestine and helping to delineate the margin around the lesion. This premixed deep sky-blue viscous solution is designed to assist with removal of polyps, adenomas and other lesions.

- Comes in prefilled syringes that are ready to use
- Outlasts saline and creates a reliable lift that will save you time

BLUEBOOST STERILE

Order Number	REF/Catalog Number	Description	Package Units
SL62521	SML10	BlueBoost Submucosal Injectable Solution, 10 ml	12/Box
SL62511	SML05	BlueBoost Submucosal Injectable Solution, 5 ml	12/Box

ELECTROSURGICAL PRECISION
GOLDKNIFE

The Micro-Tech Electrosurgical Knives feature a coated surface on the exposed metal posts that may help reduce tissue residue accumulation and provide a smooth, durable cutting surface.

The single-use T-tip GoldKnife™ has integrated electrosurgical and fluid injection functions, eliminating the need for instrument exchange which may save you time.







GOLDKNIFE T STERILE

UPN/Order Number	REF/Catalog Number	Description	Working Channel (mm)	Length (cm)	Cutting knife length when extended (mm)	Package Units
MK55131	MK-T-1-195	GoldKnife T, Irrigating	2.8	195	1.5	1/Box
MK55141	MK-T-1-235	GoldKnife T, Irrigating	2.8	235	1.5	1/Box
MK55151	MK-T-2-195	GoldKnife T, Irrigating	2.8	195	2.0	1/Box
MK55161	MK-T-2-235	GoldKnife T, Irrigating	2.8	235	2.0	1/Box

ACTIVE CORD

UPN/Order Number	REF/Catalog Number	Description	Package Units
QB51011	AEC-03-08	3.0 mm Active Cord, Reusable, Bovie 8.0 mm Male Fitting	1/Box

ELECTROSURGICAL PRECISION

GOLDKNIFE IT

The single-use Insulated Tip Injectable GoldKnife features a ceramic tip with a lumen for irrigation, that provides insulation at the distal end of the knife, minimizing invasiveness while providing a 4.0 mm cutting electrode and irrigation capability.





GOLDKNIFE IT STERILE

UPN/Order Number	REF/Catalog Number	Description	Working Channel (mm)	Working Length (cm)	Cutting Length When Extended (mm)	Package Units
MK55211	MK-IT-4-195-N	GoldKnife IT	2.8	195	4.0	1/Box
MK55221	MK-IT-4-235-N	GoldKnife IT	2.8	235	4.0	1/Box

ACTIVE CORD

UPN/Order Number	REF/Catalog Number	Description	Package Units
QB51011	AEC-03-08	3.0 mm Active Cord, Reusable, Bovie 8.0 mm Male Fitting	1/Box


TARGETED THERMAL HEMOSTASIS

ENSURE COAGULATION FORCEPS




EnSure™ reliable hemostasis of bleeding vessels without affecting surrounding tissue.

The EnSure Coagulation Forceps are designed to cauterize and coagulate GI bleeds. These single-use hemostatic forceps deliver targeted monopolar coagulation to isolate vessels and achieve hemostasis. Bleeding vessels can be approached with precision due to the reliable 360° bidirectional rotation.



Markings on catheter to help determine insertion length within the endoscope.



360° bidirectional rotation for targeting bleeds.

ENSURE COAGULATION FORCEPS STERILE

UPN/Order Number	REF/Catalog Number	Description	Jaw Style	Working Length (cm)	Outer Diameter (mm)	Jaw Opening (mm)	Package Units
CF48331	CF165-B	Ensure Coag Forceps Grooved Jaw	Grooved Jaw	165	2.3	6.3	10/Box
CF48361	CF230-B	Ensure Coag Forceps Grooved Jaw	Grooved Jaw	230	2.3	6.3	10/Box

ACTIVE CORD

UPN/Order Number	REF/Catalog Number	Description	Package Units
QB51011	AEC-03-08	3.0 mm Active Cord, Reusable, Bovie 8.0 mm Male Fitting	1/Box

GRASP, APPROXIMATE, AND CLOSE

DUAL ACTION TISSUE CLOSURE DEVICE



Gastrointestinal defect closure devices should be quick, easy to use, and effective.

The DAT™ Closure Device has a center post flanked by two arms which can operate independently to capture tissue and support manipulation to close mucosal defects.


The DAT closure device can be repositioned before deployment and can be used alone or to approximate tissue margins and hold them in place while other through-the-scope clips are placed to further secure defect closure.

Now indicated for up to 5 cm!



DAT CLOSURE DEVICE STERILE

UPN/Order Number	REF/Catalog Number	Description	Opening Width (mm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
TC35531	STA00003	DAT Closure Device	15	235	3.2	Each



DAT SLIM CLOSURE DEVICE

CLOSE MORE. REACH FURTHER.

Large defect closure demands precision, control, and access. The DAT™ Slim Closure Device builds on the proven dual-action design physicians rely on, now with 360° rotation and compatibility with a standard 2.8 mm working channel.

DAT Slim can be repositioned before deployment and used alone or to approximate tissue margins while additional clips are placed to further secure defect closure.

- 360° rotation via catheter (by the physician) or handle (by the nurse/tech) for precise alignment
- Compatible with a standard 2.8 mm working channel
- Available in three sizes: Small, Medium, and Large
- Teeth on arms for atraumatic tissue grip and reduced slippage
- Unrestricted repositioning prior to deployment
- Indicated for mucosal and submucosal defects up to 5 cm




DAT SLIM CLOSURE DEVICE

UPN/Order Number	REF/Catalog Number	Description	Opening Width (mm)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
DA35731	DAT01002	DAT Slim Closure Device, Small	9-12	235	2.8	Each
DA35791	DAT01008	DAT Slim Closure Device, Medium	12-15	235	2.8	Each
DA35851	DAT0100E	DAT Slim Closure Device, Large	15-18	235	2.8	Each


LOCK INTO TISSUE

LOCKADO HEMOSTASIS CLIP

Lockado® allows unrestricted opening and closing so it can be repositioned until the desired position is achieved.



Hemostasis clips need to provide accuracy, reliability and precise functionality. With two sets of opposing atraumatic teeth in the arms of the clip, Lockado delivers a reliable and secure grip for consistent performance. Driven by Micro-Tech Endoscopy's proven clip delivery system, the Lockado clip provides exceptional reliability and tissue management.






LOCKADO MAX HEMOSTASIS CLIP* STERILE

UPN/Order Number	REF/Catalog Number	Opening Width (mm)	Sheath Diameter (mm)	Working Length (cm)	Minimum Channel Size (mm)	Packing Units
LO35041	LOCK-D-26-235-C-N	22	2.6	235	2.8	10/Box
LO35045	LOCK-D-26-235-C-N	22	2.6	235	2.8	2/Box

*Jaws may not open to maximum reach until placed in contact with tissue.




PROTECT TISSUE SAMPLES
ESD SPECIMEN TRAY

This convenient specimen container has everything needed to mount tissue specimens before sending them to pathology. The tabs on the lid help seal the box so tissue is protected while it is transported. A tube of push pins is included in the box for mounting the tissue to the specimen pad. Markings on the blue pad help track orientation and include a measurement grid for sizing.

ESD SPECIMEN TRAY

UPN/Order Number	REF/Catalog Number	Description	Dimensions	Package Units
QB55011	ESDT-9090	ESD Specimen Tray with pins	5.5″ x 5.5″ x 2.75″	6/Box



ORGANIZE, PROTECT AND SECURE
DEVICE ORGANIZER

The Micro-Tech device organizer is perfect for storing medical devices during advanced endoscopic procedures. Without organization, devices may fall to the floor, get tangled or even become damaged. The organizer conveniently adheres to surfaces such as the countertop, bedside rail, or endoscopy cart.

Keep devices organized, protected and secure until needed.

- Adhesive anchoring strips
- Five oversized, clear pouches
- Bendable hooks to hold ring handles

MICRO-TECH DEVICE ORGANIZER

UPN/Order Number	REF/Catalog Number	Description	Dimensions (cm)	Package Units
QL22951	EDO-1000-5	Endoscopy Device Organizer, 5 Pocket	100 x 100	20/Box

POLYPECTOMY-RETRIEVAL

POLYPECTOMY

BRAIDED POLYPECTOMY SNARES

FIRM POLYPECTOMY SNARES

COLD POLYPECTOMY SNARES

ENDOINK

BLUE BEACON

RETRIEVAL

RETRIEVAL FORCEPS


PROTRAP

POLYP TRAPS

ENTRAP

PRECISE CUT

POLYPECTOMY SNARES



Our line of hot and cold polypectomy snares features effective cutting wires and user-friendly handles at cost-effective prices

- With the innovative power connector, the snare can be operated and rotated with ease
- Our hot snares are available in a compliant braided wire or a firmer wire for effective tissue acquisition and resection



KEY BENEFITS

SHORTER THROW HANDLE

The handle throw is adjusted according to the snare size. Small and Medium snares have a shorter throw than Large or Extra Large snares. A shorter throw improves comfort for users with smaller hands.

BRAIDED WIRE

The braided wire combines the tissue acquisition features of a firm snare with the flexibility and smooth deployment of a soft snare. Suitable for daily use in facilities looking to minimize inventory.

INNOVATIVE POWER CONNECTION

The handle features a dedicated power connection, separate from the finger ring slider. This feature permits unrestricted movement and rotation of the snare handle and may reduce procedural delays.

BRAIDED SOFT SNARES STERILE

UPN/Order Number	REF/Catalog Number	Description	Loop Diameter (mm)	Sheath Diameter (Fr)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
PS51011	NPFS01-01023230	Mini Braided Oval	10	7	230	2.8	10/Box
PS51021	NPFS01-01523230	Small Braided Oval	15	7	230	2.8	10/Box
PS51031	NPFS01-02023230	Medium Braided Oval	20	7	230	2.8	10/Box
PS51041	NPFS01-02423230	Large Braided Oval	24	7	230	2.8	10/Box
PS51051	NPFS01-03023230	Extra Large Braided Oval	30	7	230	2.8	10/Box

FIRM SNARES STERILE

UPN/Order Number	REF/Catalog Number	Description	Loop Diameter (mm)	Sheath Diameter (Fr)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
PS51071	NPFS02-01023230	Mini Oval Firm	10	7	230	2.8	10/Box
PS51081	NPFS02-01523230	Small Oval Firm	15	7	230	2.8	10/Box
PS51091	NPFS02-02023230	Medium Oval Firm	20	7	230	2.8	10/Box
PS51101	NPFS02-02423230	Large Oval Firm	24	7	230	2.8	10/Box
PS51111	NPFS02-03023230	Extra Large Oval Firm	30	7	230	2.8	10/Box
PS51121	NPFS02-03623230	Jumbo Oval Firm	36	7	230	2.8	10/Box

FIRM SNARES HEX STERILE

UPN/Order Number	REF/Catalog Number	Description	Loop Diameter (mm)	Sheath Diameter (Fr)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
PS52031	NPFS02-11523230	Small Hex Firm	15	7	230	2.8	10/Box
PS52041	NPFS02-12523230	Large Hex Firm	25	7	230	2.8	10/Box

ACTIVE CORD

UPN/Order Number	REF/Catalog Number	Description	Package Units
QB51011	AEC-03-08	3.0 mm Active Cord, Reusable, Bovie 8.0 mm Male Fitting	1/Box



EFFICIENT RESECTION

COLD SNARES

- Designed specifically to resect diminutive polyps
- Available in two sizes – 10 mm and 15 mm
- Thin cutting wire provides for a guillotine-type cut
- Loop shape assists in capturing smaller polyps
- Stiffer catheter aids in insertion and provides support during resection



KEY BENEFITS

DEDICATED COLD SNARE

Specifically designed for resecting diminutive polyps using a cold snare technique, the oval or diamond shape loop and flexible wire, capture smaller polyps with ease.

THIN WIRE

The thin .23 mm cutting wire provides a guillotine cut that reduces the force needed to complete the resection and may limit tearing or snagging.

COMFORTABLE ERGONOMIC HANDLES

Accommodating various hand sizes, the ergonomic three-ring handle features a shortened throw and comfortable grip. The stiff catheter provides support for a smooth effortless cut.

COLD SNARES STERILE

UPN/Order Number	REF/Catalog Number	Description	Loop Diameter (mm)	Sheath Diameter (Fr)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
CS50011	CS3-11023230	Diamond Thin Wire	10	7	230	2.8	10/Box
CS50021	CS3-11523230	Diamond Thin Wire	15	7	230	2.8	10/Box
CS50031	CS3-21023230	Oval Thin Wire	10	7	230	2.8	10/Box
CS50041	CS3-21523230	Oval Thin Wire	15	7	230	2.8	10/Box

ENDOSCOPIC TATTOO
ENDOINK

EndoINK™ Endoscopic Marker is a single-use suspension fluid designed to permanently mark lesions in the GI tract. Developed from high purity carbon particles, EndoINK remains in the submucosal injection site, creating a permanent mark to facilitate surgical removal and examination.



KEY BENEFITS

RELIABLE MARKER
Endoscopic tattoo for reliable marking of lesions and areas of reference in the GI tract.

SYRINGE CONFIGURATIONS
Conveniently available in 5 ml and 10 ml configuration to accommodate catheter priming.

EASY TO USE
Sterile fluid supplied in a capped syringe for single patient use.

ENDOINK STERILE ②


UPN/Order Number	REF/Catalog Number	Description	Package Units
NK62711	ENK-05-05	EndoINK Endoscopic Marker, 5 ml	5/Box
NK62715	ENK-05-12	EndoINK Endoscopic Marker, 5 ml	12/Box
NK62721	ENK-10-05	EndoINK Endoscopic Marker, 10 ml	5/Box
NK62725	ENK-10-12	EndoINK Endoscopic Marker, 10 ml	12/Box

TRUE BLUE

BLUE BEACON GI MARKER



- Blue Beacon™ is useful in GI procedures that require a contrast stain to temporarily mark or highlight tissue
- Can be diluted and introduced to the irrigation, sprayed or injected around the area of interest
- The 0.4% Indigo Carmine solution is supplied sterile in a convenient luer lock 5 ml syringe, 5 syringes per box
- The 0.8% Indigo Carmine solution is supplied sterile in a 10 ml single-use vial, 10 vials per box



BLUE BEACON GI MARKER STERILE

UPN/Order Number	REF/Catalog Number	Description	Package Units
IC62021	P-DYZ-04-05	Blue Beacon GI Marker 5 ml Syringe Indigo Carmine 0.4% Sterile	5/Box
IC62061	DYZ-08-10	Blue Beacon GI Marker 10 ml Vial Indigo Carmine 0.8% Sterile	10/Box

This solution is NOT FOR INTRAVENOUS OR INTRAMUSCULAR INJECTION.
It is a single-use device indicated for temporary marking of the GI tract.

SPRAY CATHETER STERILE

UPN/Order Number	REF/Catalog Number	Description	Working Length (cm)	Minimum Channel Size (mm)	Package Units
SC62241	WP-18/2200	Spray Catheter	220	2.2	20/Box
SC62301	WP-23/2500	Spray Catheter	250	2.8	20/Box

PRECISION & CONTROL

RETRIEVAL FORCEPS



Single-use grasping forceps for secure removal of tissue or foreign bodies.

- 3 jaw types
- Ergonomic 3-ring handle
- Entry-exit markers
- Sterile single-use convenience

RETRIEVAL FORCEPS STERILE

UPN/Order Number	REF/Catalog Number	Description	Jaw Opening (mm)	Sheath Color	Working Length (cm)	Minimum Channel Size (mm)	Package Units
GF61021	NGF03-11023230	Rat Tooth	6.3	Blue	230	2.8	10/Box
GF61041	NGF13-11023230	Alligator Tooth	8.1	Blue	230	2.8	10/Box
GF61061	NGF33-11023230	Mixed Tooth	8.1	Blue	230	2.8	10/Box



ILLUMINATE | MAGNIFY | IDENTIFY

PROTRAP LUMA

Low-light environments can make identifying diminutive polyps a difficult task. The ProTrap Luma™ In-Line Polyp Trap comes with a built-in light source and 2x magnification to help visualization.

The quick-change trays allow for fast and easy tray swaps to minimize suction downtime. The trays feature a bi-directional design, which allows them to be flipped over to assist in debris clearing. When the sample is collected, the tray with sample can be directly placed in the specimen jar to minimize risks of tissue damage and biohazard exposure.

· Identify samples with confidence

· Reduce time off suction

· Reduce sample handling

PROTRAP LUMA

UPN/Order Number	REF/Catalog Number	Description	Package Units
TP54411	BBPA00001	ProTrap Luma In-Line Polyp Trap	25/Box
TP54421	BBPA10001	ProTrap Luma In-Line Polyp Trap Specimen Tray	50/Box



CLEAN COLLECTION

POLYP RETRIEVAL TRAPS

Polyp Retrieval Traps offer improved efficiency, reduced handling and cost savings. The unique Quad trap features removable collection trays and a secondary filter for peace of mind. The In-line suction polyp trap provides simple in-line installation and easy removal of polyp specimens.

POLYP TRAPS


Polyp Retrieval Traps offer improved efficiency, reduced handling and cost savings.
The unique Quad trap features removable collection trays and a secondary filter for peace of mind.

UPN/Order Number	REF/Catalog Number	Description	Package Units
QL54012	MDPT1A0002A	In-Line Suction Polyp Trap	50/Box
QL54021	PT50-50-S	Polyp Trap Single Chamber, Bulk	50/Box
PC54031	BBP-50-M	Quad Polyp Trap with Removable Chambers	25/Box

DURABLE AND GENTLE

ENTRAP POLYP RETRIEVAL NET

The Entrap™ retrieval system delivers a quick and efficient retrieval solution for polyps and tissue for histopathologic review. Utilizing advanced materials and production techniques, the Entrap facilitates improved consistency and control of the net deployment, position, and closure.



KEY BENEFITS

REFINED NET
Improved visibility through the net material for identification of the retrieval target. Durable enough for repeated deployment, but soft to gently capture and preserve the specimen.

ADVANCED PROCESSES
Cinching, tearing, and binding of the net material is reduced by the minimal net selvage and bonded net attachment. This increases efficiency and minimizes the time required to reliably capture specimens.

SMOOTH AND STRONG
A stiffer catheter material improves device insertion into the endoscope and positioning of the net. The additional columnar strength allows for gentle movement of a specimen for optimal capture.

ENTRAP RETRIEVAL SYSTEM

UPN/Order Number	REF/Catalog Number	Description	Loop Diameter (mm)	Sheath Diameter (Fr)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
NT54511	NT-E-30/23-P	Entrap Polyp Retrieval Net	30	7	230	2.8	10/Box
NT54515	NT-E-30/23-P	Entrap Polyp Retrieval Net	30	7	230	2.8	2/Box

NETIS FB RETRIEVAL NETS STERILE

UPN/Order Number	REF/Catalog Number	Description	Loop Diameter (mm)	Sheath Diameter (Fr)	Working Length (cm)	Minimum Channel Size (mm)	Package Units
MED-195-NET	MED-195-NET	Mini FB Retrieval Net	25	7.5	230	2.8	10/Box
MED-201-NET	MED-201-NET	ESD Retrieval Net	50	7.5	230	2.8	10/Box
MED-203-NET	MED-203-NET	Std FB Rotatable Net	30	7.5	230	2.8	10/Box
MED-204-NET	MED-204-NET	Jumbo Rotatable Net	40	7.5	230	2.8	10/Box
MED-217-NET	MED-217-NET	Entero FB Retrieval Net	30	7.5	260	2.8	10/Box

PRODUCT REFERENCE CHART

Order Number	Ref #	GTIN	UPN Pouch	Product Family
GW10011	MTN-BM-63/26-A	16932503510018	GW10010	Biliary
GW10021	MTN-BM-63/26-A-J	16932503510025	GW10020	Biliary
GW10031	MTN-BM-63/45-A	16932503510032	GW10030	Biliary
GW10041	MTN-BM-63/45-A-J	16932503510049	GW10040	Biliary
GW10051	MTN-BM-89/26-A	16932503510056	GW10050	Biliary
GW10071	MTN-BM-89/26-A-J	16932503510070	GW10070	Biliary
GW10091	MTN-BM-89/45-A	16932503510094	GW10090	Biliary
GW10111	MTN-BM-89/45-A-J	16932503510117	GW10110	Biliary
GW10061	MTN-BM-89/26-A-S	16932503510063	GW10060	Biliary
GW10081	MTN-BM-89/26-A-J-S	16932503510087	GW10080	Biliary
GW10101	MTN-BM-89/45-A-S	16932503510100	GW10100	Biliary
GW10121	MTN-BM-89/45-A-J-S	16932503510124	GW10120	Biliary
SS75761	RDSP-0520-221-1	16932503575765	SS75760	Biliary
SS75781	RDSP-0520-231-1	16932503575789	SS75780	Biliary
SS16121	RDSP-0525-221-1	16932503516126	SS16120	Biliary
SS16141	RDSP-0525-231-1	16932503516140	SS16140	Biliary
SS16181	RDSP-0530-221-1	16932503516188	SS16180	Biliary
SS16201	RDSP-0530-231-1	16932503516201	SS16200	Biliary
SS75041	RDSP-0520-221	16932503575048	SS75040	Biliary
SS75061	RDSP-0520-231	16932503575062	SS75060	Biliary
SS75101	RDSP-0525-221	16932503575109	SS75100	Biliary
SS75121	RDSP-0525-231	16932503575123	SS75120	Biliary
SS75161	RDSP-0530-221	16932503575161	SS75160	Biliary
SS75181	RDSP-0530-231	16932503575185	SS75180	Biliary
SS75741	RDSP-0520-211-1	16932503575741	SS75740	Biliary
SS75801	RDSP-0525-211-1	16932503575802	SS75800	Biliary
SS16161	RDSP-0530-211-1	16932503516164	SS16160	Biliary
SS75021	RDSP-0520-211	16932503575024	SS75020	Biliary
SS75081	RDSP-0525-211	06932503575089	SS75080	Biliary
SS75141	RDSP-0530-211	06932503575140	SS75140	Biliary
SB75921	RSRB-P-9/12/15-20	06932503575928	SB75920	Biliary
SB75931	RSRB-P-12/15/18-20	06932503575935	SB75930	Biliary
SB75941	RSRB-P-15/18/21-20	16932503575949	SB75940	Biliary
SB76131	RSRB-P-9/12/15-20-B	16932503576137	SB76130	Biliary
SB76141	RSRB-P-12/15/18-20-B	16932503576144	SB76140	Biliary
SB76151	RSRB-P-15/18/21-20-B	16932503576151	SB76150	Biliary
GL62781	MT-RGL-O	16932503562789	GL62780	Biliary
GL62791	MT-RGL-P	16932503562796	GL62790	Biliary
BF40661	SCNBF13-10022180	16932503540664	BF40660	Biopsy
BF40641	SCNBF13-10122180	16932503540640	BF40640	Biopsy
BF40671	SCNBF13-10022230	16932503540671	BF40670	Biopsy
BF40676	SCNBF13-10022230	66932503540676	BF40674	Biopsy
BF40651	SCNBF13-10122230	16932503540657	BF40650	Biopsy
BF40656	SCNBF13-10122230	66932503540652	BF40654	Biopsy
BF40271	LCNBF13-10024180	16932503540275	BF40270	Biopsy
BF40291	LCNBF13-10124180	16932503540299	BF40290	Biopsy
BF40281	LCNBF13-10024230	16932503540282	BF40280	Biopsy
BF40286	LCNBF13-10024230	66932503540287	BF40284	Biopsy
BF40301	LCNBF13-10124230	16932503540305	BF40300	Biopsy
BF40306	LCNBF13-10124230	66932503540300	BF40304	Biopsy

Order Number	Ref #	GTIN	UPN Pouch	Product Family
BF40461	JSNBF63-10028230	16932503540466	BF40460	Biopsy
BF40451	JSNBF63-10128230	16932503540459	BF40450	Biopsy
BF40491	JSNBF63-10133230	16932503540497	BF40490	Biopsy
BF40021	NBF02-11018160	16932503540022	BF40020	Biopsy
BF40051	NBF02-11118160	16932503540053	BF40050	Biopsy
BF40141	NBF12-11018160	16932503540145	BF40140	Biopsy
BF40401	NBF12-11118160	16932503540404	BF40400	Biopsy
BF40091	NBF03-11023230	16932503540091	BF40090	Biopsy
BF40121	NBF03-11123230	16932503540121	BF40120	Biopsy
BF40096	NBF03-11023230	66932503540096	BF40094	Biopsy
BF40126	NBF03-11123230	66932503540126	BF40124	Biopsy
BF40201	NBF13-11023230	16932503540206	BF40200	Biopsy
BF40231	NBF13-11123230	16932503540237	BF40230	Biopsy
BF40206	NBF13-11023230	66932503540201	BF40200	Biopsy
BF40236	NBF13-11123230	66932503540232	BF40234	Biopsy
BF40441	NBF63-11023230	16932503540442	BF40440	Biopsy
BF40421	NBF63-11123230	16932503540428	BF40420	Biopsy
BF40446	NBF63-11023230	66932503540447	BF40444	Biopsy
BF40426	NBF63-11123230	66932503540423	BF40424	Biopsy
CY49021	CYB-20S-18C	16932503549025	CY49020	Biopsy
CY49051	CYB-30S-24C	16932503549056	CY49050	Biopsy
MC924A-10	MC924A-10	10666416140333	MC924A	Biopsy
BF51811	000381	16942852451818	BF51810	Biopsy
BF51821	000382	16942852451825	BF51820	Biopsy
BF51831	000383	16942852451832	BF51830	Biopsy
BF51841	000384	16942852451849	BF51840	Biopsy
BF51851	000385	16942852451856	BF51850	Biopsy
BF51861	000386	16942852451863	BF51860	Biopsy
BF51881	000389	16942852451887	BF51880	Biopsy
BF51871	000388	16942852451870	BF51870	Biopsy
BF63851	00389B	16942852463835	BF63830	Biopsy
BF63831	00386B	16942852463859	BF63850	Biopsy
BF52071	00821C	16942852452075	BF52070	Biopsy
BF52091	00823C	16942852452099	BF52090	Biopsy
BF52141	100512	16942852452143	BF52140	Biopsy
BF52151	100513	16942852452150	BF52150	Biopsy
CD63681	CDS11007	16932503563687	CD63680	Cholangioscopy
CD63701	CDS11008	16932503563700	CD63700	Cholangioscopy
CD63691	CDS22007	16932503563694	CD63690	Cholangioscopy
CD63761	CDS22008	16932503563762	CD63760	Cholangioscopy
SW71212	BS-W-100	26932503571214	SW71210	Cholangioscopy
DF49425	BF10006	56932503549429	DF49420	Cholangioscopy
DF49365	BF16006	56932503549368	DF49360	Cholangioscopy
BE11721	CEB00000	16932503511725	BE11720	Cholangioscopy
BE11731	CEB01010	16932503511732	BE11730	Cholangioscopy
CD63241	CDS22005	16932503563243	CD63240	Cholangioscopy
CD63251	CDS22005	16932503563250	CD63250	Cholangioscopy
BE11861	CDS22005	16932503511862	BE11860	Cholangioscopy
BE11871	CDS22005	16932503511879	BE11870	Cholangioscopy
AM71601	CDS33101	16932503571606	AM71600	Cholangioscopy
AM71611	CDS44101	16932503571613	AM71610	Cholangioscopy
AM71631	CDS66101	16932503571637	AM71630	Cholangioscopy
EL27C115VEHL	EL27C115VEHL	04260667540011	EL27C115VEHL	Cholangioscopy

Order Number	Ref #	GTIN	UPN Pouch	Product Family
P243000FST	P243000FST	04260667540097	P243000FST	Cholangioscopy
P303000FST	P303000FST	04260667540158	P303000FST	Cholangioscopy
P453000FST	P453000FST	04260667540196	P453000FST	Cholangioscopy
MB23031	MBD-AM-BB-2	16932503523032	MB23030	Dilation
MB23071	MBD-BM-BB-2	16932503523070	MB23070	Dilation
MB23111	MBD-CM-BB-2	16932503523117	MB23110	Dilation
MB23151	MBD-DM-BB-2	16932503523155	MB23150	Dilation
MB23191	MBD-EM-BB-2	16932503523193	MB23190	Dilation
MB23221	MBD-FM-BB-2	16932503523223	MB23220	Dilation
FW24015	FWB-06/80-6/18-M	56932503524013	FW24014	Dilation
FW24025	FWB-08/80-6/18-M	56932503524020	FW24024	Dilation
FW24035	FWB-10/80-6/18-M	56932503524037	FW24034	Dilation
FW24045	FWB-12/80-6/18-M	56932503524044	FW24044	Dilation
FW24055	FWB-15/80-6/18-M	56932503524051	FW24054	Dilation
FW24065	FWB-18/80-6/18-M	56932503524068	FW24064	Dilation
FW24165	FWB-06/55-6/23-S	56932503524167	FW24164	Dilation
FW24175	FWB-08/55-6/23-S	56932503524174	FW24174	Dilation
FW24185	FWB-10/55-6/23-S	56932503524181	FW24184	Dilation
FW24195	FWB-12/55-6/23-S	56932503524198	FW24194	Dilation
FW24205	FWB-15/55-6/23-S	56932503524204	FW24204	Dilation
FW24215	FWB-18/55-6/23-S	56932503524211	FW24214	Dilation
FW24225	FWB-20/55-6/23-S	56932503524228	FW24224	Dilation
FW24075	FWB-06/80-6/18-S	56932503524075	FW24074	Dilation
FW24085	FWB-08/80-6/18-S	56932503524082	FW24084	Dilation
FW24095	FWB-10/80-6/18-S	56932503524099	FW24094	Dilation
FW24105	FWB-12/80-6/18-S	56932503524105	FW24104	Dilation
FW24125	FWB-15/80-6/18-S	56932503524129	FW24124	Dilation
FW24145	FWB-18/80-6/18-S	56932503524143	FW24144	Dilation
FW24155	FWB-20/80-6/18-S	56932503524150	FW24154	Dilation
DB25011	BID-15-60-01	16932503525012	DB25010	Dilation
DB25012	BID-15-60-01	26932503525019		Dilation
DG25111	DG-18-08-21	16932503525111	DG25110	Dilation
CB62411	CB10-0030/09-S	16932503562413	CB62410	Dilation
ES34621	EUS-25-0-N	16932503534625	ES34620	EUS
ES34601	EUS-22-0-N	16932503534601	ES34600	EUS
ES34641	EUS-22-0-F	16932503534649	ES34640	EUS
ES34631	EUS-25-1-N	16932503534632	ES34630	EUS
ES34611	EUS-25-1-N	16932503534618	ES34610	EUS
ES34591	EUS-22-1-N	16932503534595	ES34590	EUS
ES34671	EUS-22-1-F	16932503534670	ES34670	EUS
ES34661	EUS-22-1-F	16932503534663	ES34660	EUS
LO30861	LOCK-D-26-235-C-N	16932503530863	LO30860	Hemostasis
LO30865	LOCK-D-26-235-C-N	56932503530861	LO30864	Hemostasis
LO30891	LOCK-F-26-235-C-N	16932503530894	LO30890	Hemostasis
LO30895	LOCK-F-26-235-C-N	56932503530892	LO30894	Hemostasis
LO35041	LOCK-G-26-235-C-N	16932503535042	LO35040	Hemostasis
LO35045	LOCK-G-26-235-C-N	56932503535040	LO35044	Hemostasis
RC40011	DC0235	26932619740016	RC40010	Hemostasis
RC40031	DC0235W	26932619740030	RC40030	Hemostasis
RC30415	ROCC-C-26-235-C-N	56932503530410	RC30414	Hemostasis
RC30411	ROCC-C-26-235-C-N	16932503530412	RC30410	Hemostasis
RC30445	ROCC-D-26-235-C-N	56932503530441	RC30444	Hemostasis
RC30441	ROCC-D-26-235-C-N	16932503530443	RC30440	Hemostasis
RC30385	ROCC-F-26-235-C-R	56932503530380	RC30384	Hemostasis

Order Number	Ref #	GTIN	UPN Pouch	Product Family
RC30381	ROCC-F-26-235-C-R	16932503530382	RC30380	Hemostasis
RC30545	ROCC-H-26-235-C-N	56932503530540	RC30544	Hemostasis
RC30541	ROCC-H-26-235-C-N	16932503530542	RC30540	Hemostasis
RT30665	ROCC-D-26-235-C	56932503530663	RT30664	Hemostasis
RT30661	ROCC-D-26-235-C	16932503530665	RT30660	Hemostasis
RT30785	ROCC-F-26-235-C	56932503530786	RT30784	Hemostasis
RT30781	ROCC-F-26-235-C	16932503530788	RT30780	Hemostasis
IN33301	IN12-223232302	16932503533307	IN33300	Hemostasis
IN33331	IN12-223232302	16932503533338	IN33330	Hemostasis
IN33131	IN12-225231802	16932503533130	IN33130	Hemostasis
IN33221	IN12-225231802	16932503533222	IN33220	Hemostasis
IN33151	IN12-225232302	16932503533154	IN33150	Hemostasis
IN33241	IN12-225232302	16932503533246	IN33240	Hemostasis
ML32931	MBLS-6-NL-S	16932503532935	ML32930	Hemostasis
HL32161	HMBLS-4F-NL	16932503532164	HL32160	Hemostasis
ML32651	MBLS-6F-NL	16932503532652	ML32650	Hemostasis
QL30011	MDHL100000A	00858276007410	QL30010	Hemostasis
QL71151	AW01-OLY3	00850025860598	QL71150	Infection Prevention
QL71101	AW01-OLY4	00850025860376	QL71100	Infection Prevention
QL71112	AW01-OLY5	00850025860321		Infection Prevention
QL70021	MD-A-BV02	00858276007786	QL70020	Infection Prevention
QL70051	MD-B-BV02	00858276007359	QL70050	Infection Prevention
QL72191	BW01-OLYB	10850025860748	QL72190	Infection Prevention
QL72201	BW01-PENB	10850025860762	QL72200	Infection Prevention
QL72111	BW01-OLY5	00850025860550	QL72110	Infection Prevention
QL72101	BW01-OLY4	10850025860724	QL72100	Infection Prevention
QL72151	BW01-OLY3	10850025860731	QL72150	Infection Prevention
QL72191	BW01-OLYB	10850025860748	QL72190	Infection Prevention
QL72171	BW01-OLYA	10850025860656	QL72170	Infection Prevention
QL72231	BW01-OLYC	10850025860663	QL72230	Infection Prevention
QL72131	BW01-FUJ5-F7	00850025860512	QL72130	Infection Prevention
QL72281	BW01-FUJ4-F7	10850025860717	QL72280	Infection Prevention
QL72301	BW01-FUJ3-F7	10850025860878	QL72300	Infection Prevention
QL72211	BW01-FUJA	10850025860649	QL72210	Infection Prevention
QL72261	BW01-FUJC-F7	10850025860694	QL72260	Infection Prevention
QL72121	BW01-PEN4	00850025860574	QL72120	Infection Prevention
QL72161	BW01-PEN4	10850025860755	QL72160	Infection Prevention
QL72201	BW01-PENB	10850025860762	QL72200	Infection Prevention
QL72181	BW01-PENA	10850025860670	QL72180	Infection Prevention
QL76101	BW02-OLY	10850025860779	QL76100	Infection Prevention
QL76111	BW02-OLYC02	00850025860789	QL76110	Infection Prevention
QL76121	BW02-IRRTUB	10850025860793	QL76120	Infection Prevention
QL76151	BW02-IRRTUBSL	10850025860847	QL76150	Infection Prevention
QL76311	BW02-CO2TUBOLY	10858276007974	QL76310	Infection Prevention
QL76321	BW02-CO2TUBAOLY	10850025860908	QL76320	Infection Prevention
QL76211	BW02-CO2CONOLY	10850025860854	QL76210	Infection Prevention
QL76221	BW02-CO2CONFJ	10850025860861	QL76220	Infection Prevention
QL26031	MDBB1AY003A	00850025860048	QL26030	Infection Prevention
BB60011	AC01-103.A	16932503560013	BB60010	Infection Prevention
BB60021	AC01-103.P	16932503560020	BB60020	Infection Prevention
BB60041	AC01-108.M	16932503560044	BB60040	Infection Prevention
BB60051	AC01-109.M	16932503560051	BB60050	Infection Prevention
QL78111	BW03-DTPL	10850025860816	QL78110	Infection Prevention
QL78101	BW03-DTPS	10850025860809	QL78100	Infection Prevention

Order Number	Ref #	GTIN	UPN Pouch	Product Family
CP62131	CTP-02025	16932503562130	CP62130	Infection Prevention
QL22011	MDBR1A0504A	00858276007014	QL22010	Infection Prevention
QL22021	MDBR1B0016A	00858276007045	QL22020	Infection Prevention
QL22031	MDBR1B0015A	00858276007076	QL22030	Infection Prevention
QL22041	MDBR1B0005A	00858276007106	QL22040	Infection Prevention
QL22051	MDBR1B0013A	00858276007137	QL22050	Infection Prevention
QL22061	MDBR1B0014A	00858276007168	QL22060	Infection Prevention
QL22561	KT1A0006A	00850025860222	QL22560	Infection Prevention
QL22511	MDKT1A0001A	00858276007441	QL22510	Infection Prevention
QL22551	MDKT1A0005A	00850025860079	QL22550	Infection Prevention
QL22531	MDKT1A0003A	00858276007502	QL22530	Infection Prevention
QL22711	MDJL1A0001A	00858276007502	QL22530	Infection Prevention
QL22811	MDSB1B0001A	00858276007656	QL22810	Infection Prevention
ET80010	NST01-222-18.060	06932503580014	ET80010	Nitinol Stents
ET80020	NST01-222-18.080	06932503580021	ET80020	Nitinol Stents
ET80030	NST01-222-18.100	06932503580038	ET80030	Nitinol Stents
ET80040	NST01-222-18.120	06932503580045	ET80040	Nitinol Stents
ET80050	NST01-222-18.140	06932503580052	ET80050	Nitinol Stents
ET80060	NST01-222-18.150	06932503580069	ET80060	Nitinol Stents
ET80070	NST01-222-20.060	06932503580076	ET80070	Nitinol Stents
ET80080	NST01-222-20.080	06932503580083	ET80080	Nitinol Stents
ET80090	NST01-222-20.100	06932503580090	ET80090	Nitinol Stents
ET80100	NST01-222-20.120	06932503580106	ET80100	Nitinol Stents
ET80110	NST01-222-20.140	06932503580113	ET80110	Nitinol Stents
ET80120	NST01-222-20.150	06932503580120	ET80120	Nitinol Stents
ET80130	NST01-222-22.060	06932503580137	ET80130	Nitinol Stents
ET80140	NST01-222-22.080	06932503580144	ET80140	Nitinol Stents
ET80150	NST01-222-22.100	06932503580151	ET80150	Nitinol Stents
ET80160	NST01-222-22.120	06932503580168	ET80160	Nitinol Stents
ET80170	NST01-222-22.140	06932503580175	ET80170	Nitinol Stents
ET80180	NST01-222-22.150	06932503580182	ET80180	Nitinol Stents
ET80190	NST01-224-18.060	06932503580199	ET80190	Nitinol Stents
ET80200	NST01-224-18.080	06932503580205	ET80200	Nitinol Stents
ET80210	NST01-224-18.100	06932503580212	ET80210	Nitinol Stents
ET80220	NST01-224-18.120	06932503580229	ET80220	Nitinol Stents
ET80230	NST01-224-18.140	06932503580236	ET80230	Nitinol Stents
ET80240	NST01-224-18.150	06932503580243	ET80240	Nitinol Stents
ET80250	NST01-224-20.060	06932503580250	ET80250	Nitinol Stents
ET80260	NST01-224-20.080	06932503580267	ET80260	Nitinol Stents
ET80270	NST01-224-20.100	06932503580274	ET80270	Nitinol Stents
ET80280	NST01-224-20.120	06932503580281	ET80280	Nitinol Stents
ET80290	NST01-224-20.140	06932503580298	ET80290	Nitinol Stents
ET80300	NST01-224-20.150	06932503580304	ET80300	Nitinol Stents
ET80310	NST01-224-22.060	06932503580311	ET80310	Nitinol Stents
ET80320	NST01-224-22.080	06932503580328	ET80320	Nitinol Stents
ET80330	NST01-224-22.100	06932503580335	ET80330	Nitinol Stents
ET80340	NST01-224-22.120	06932503580342	ET80340	Nitinol Stents
ET80350	NST01-224-22.140	06932503580359	ET80350	Nitinol Stents
ET80360	NST01-224-22.150	06932503580366	ET80360	Nitinol Stents
ET80510	NST61-224-18.060	06932503580519	ET80510	Nitinol Stents
ET80520	NST61-224-18.080	06932503580526	ET80520	Nitinol Stents
ET80530	NST61-224-18.100	06932503580533	ET80530	Nitinol Stents
ET80540	NST61-224-18.120	06932503580540	ET80540	Nitinol Stents

Order Number	Ref #	GTIN	UPN Pouch	Product Family
ET80550	NST61-224-18.140	06932503580557	ET80550	Nitinol Stents
ET80560	NST61-224-20.060	06932503580564	ET80560	Nitinol Stents
ET80570	NST61-224-20.080	06932503580571	ET80570	Nitinol Stents
ET80580	NST61-224-20.100	06932503580588	ET80580	Nitinol Stents
ET80590	NST61-224-20.120	06932503580595	ET80590	Nitinol Stents
ET80600	NST61-224-20.140	06932503580601	ET80600	Nitinol Stents
ET80610	NST61-224-22.060	06932503580618	ET80610	Nitinol Stents
ET80620	NST61-224-22.080	06932503580625	ET80620	Nitinol Stents
ET80630	NST61-224-22.100	06932503580632	ET80630	Nitinol Stents
ET80640	NST61-224-22.120	06932503580649	ET80640	Nitinol Stents
ET80650	NST61-224-22.140	06932503580656	ET80650	Nitinol Stents
TC35531	STA00003	16932503535530	TC35530	EMR-ESD
DA35731	DAT01002	16932503535738	DA35730	EMR-ESD
DA35791	DAT01008	16932503535790	DA35790	EMR-ESD
DA35851	DAT0100E	16932503535851	DA35850	EMR-ESD
QL22951	EDO-1000-5	00850025860192	QL22950	EMR-ESD
SL62521	SML10	16932503562529	SL62520	EMR-ESD
SL62511	SML05	16932503562512	SL62510	EMR-ESD
PS51181	PS51181	16932503551189	PS51180	EMR-ESD
PS52051	NRPFS05-01523230	16932503552056	PS52050	EMR-ESD
PS52061	NRPFS05-02023230	06932503552066	PS52060	EMR-ESD
CF48331	CF165-B	16932503548332	CF48330	EMR-ESD
CF48361	CF230-B	16932503548363	CF48360	EMR-ESD
QB51011	AEC-03-08	00850025860253	QB51010	EMR-ESD
MK55131	MK-T-1-195	16932503555132	MK55130	EMR-ESD
MK55141	MK-T-1-235	16932503555149	MK55140	EMR-ESD
MK55151	MK-T-2-195	16932503555156	MK55150	EMR-ESD
MK55161	MK-T-2-235	16932503555163	MK55160	EMR-ESD
MK55211	MK-IT-4-195-N	16932503555217	MK55210	EMR-ESD
MK55221	MK-IT-4-235-N	16932503555224	MK55220	EMR-ESD
CS50051	CS2-21023231	16932503550052	CS50050	EMR-ESD
CS50061	CS2-21523231	16932503550069	CS50060	EMR-ESD
CS50071	CS2-22023231	16932503550076	CS50070	EMR-ESD
LO35041	LOCK-D-26-235-C-N	16932503535042	LO35040	EMR-ESD
LO35045	LOCK-D-26-235-C-N	56932503535040	LO35044	EMR-ESD
TS63041	ETS00004	16932503563045		EMR-ESD
TS63061	ETS00006	16932503563069		EMR-ESD
TS63101	ETS00004	16932503563106		EMR-ESD
TS63121	ETS00006	16932503563120		EMR-ESD
PS51011	NPFS01-01023230	16932503551011	PS51010	Polypectomy & Retrieval
PS51021	NPFS01-01523230	16932503551028	PS51020	Polypectomy & Retrieval
PS51031	NPFS01-02023230	16932503551035	PS51030	Polypectomy & Retrieval
PS51041	NPFS01-02423230	16932503551042	PS51040	Polypectomy & Retrieval
PS51051	NPFS01-03023230	16932503551059	PS51050	Polypectomy & Retrieval
PS51071	NPFS02-01023230	16932503551073	PS51070	Polypectomy & Retrieval
PS51081	NPFS02-01523230	16932503551080	PS51080	Polypectomy & Retrieval
PS51091	NPFS02-02023230	16932503551097	PS51090	Polypectomy & Retrieval
PS51101	NPFS02-02423230	16932503551103	PS51100	Polypectomy & Retrieval
PS51111	NPFS02-03023230	16932503551110	PS51110	Polypectomy & Retrieval
PS51121	NPFS02-03623230	16932503551127	PS51120	Polypectomy & Retrieval
PS52031	NPFS02-11523230	16932503552032	PS52030	Polypectomy & Retrieval
PS52041	NPFS02-12523230	16932503552049	PS52040	Polypectomy & Retrieval
CS50011	CS3-11023230	16932503550014	CS50010	Polypectomy & Retrieval

Order Number	Ref #	GTIN	UPN Pouch	Product Family
CS50021	CS3-11523230	16932503550021	CS50020	Polypectomy & Retrieval
CS50031	CS3-21023230	16932503550038	CS50030	Polypectomy & Retrieval
CS50041	CS3-21523230	16932503550045	CS50040	Polypectomy & Retrieval
NK62711	ENK-05-05	16932503562710	NK62710	Polypectomy & Retrieval
NK62715	ENK-05-12	56932503562718	NK62714	Polypectomy & Retrieval
NK62721	ENK-10-05	16932503562727	NK62720	Polypectomy & Retrieval
NK62725	ENK-10-12	56932503562725	NK62724	Polypectomy & Retrieval
IC62021	P-DYZ-04-05	16932503562062	IC62060	Polypectomy & Retrieval
IC62061	DYZ-08-10	16932503562024	IC62020	Polypectomy & Retrieval
SC62241	WP-18/2200	16932503562246	SC62240	Polypectomy & Retrieval
SC62301	WP-23/2500	16932503562307	SC62300	Polypectomy & Retrieval
GF61021	NGF03-11023230	16932503561027	GF61020	Polypectomy & Retrieval
GF61041	NGF13-11023230	16932503561041	GF61040	Polypectomy & Retrieval
GF61061	NGF33-11023230	16932503561065	GF61060	Polypectomy & Retrieval
TP54411	BBPA00001	00858276007922	TP54410	Polypectomy & Retrieval
TP54421	BBPA10001	00858276007953	TP54420	Polypectomy & Retrieval
QL54012	MDPT1A0002A	00858276007199	QL54010	Polypectomy & Retrieval
QL54021	PT50-50-S	10858276007707	QL54020	Polypectomy & Retrieval
PC54031	BBP-50-M	16932503554036	PC54030	Polypectomy & Retrieval
NT54511	NT-E-30/23-P	16932503554517	NT54510	Polypectomy & Retrieval
NT54515	NT-E-30/23-P	56932503554515	NT54514	Polypectomy & Retrieval
MED-195-NET	MED-195-NET	08054377481957		Polypectomy & Retrieval
MED-201-NET	MED-201-NET	08054377483012		Polypectomy & Retrieval
MED-203-NET	MED-203-NET	08054377483258		Polypectomy & Retrieval
MED-204-NET	MED-204-NET	08054377483241		Polypectomy & Retrieval
MED-217-NET	MED-217-NET	08054377482176		Polypectomy & Retrieval

TERMS AND POLICIES

SHIPMENT AND PAYMENT TERMS

Items are normally shipped Federal Express Ground or UPS Ground, FOB shipping location. Shipping charges will be added to the invoice. Payment Terms are net 30 days. All prices are subject to change without notice. Express deliveries of stock items can be arranged and will be subject to additional freight charges.

SPECIAL ORDER ITEMS

In addition to stock items listed within this catalog, Micro-Tech Endoscopy manufactures non-stock items. Non-stock items may be subject to longer lead times. Contact Customer Support prior to ordering for more information.

STANDING ORDER SHIPMENTS

A standing order shipment simplifies ordering and shipments of regularly purchased products. Products ordered routinely on a regular basis can be scheduled to ship at an agreed upon, regular interval. Standing order shipments offer convenience and can be canceled at any time.

RETURNED GOODS POLICY

Micro-Tech Endoscopy offers a 100% refund on all unused and unopened products in the original product box, that are returned within 90 days from the date of purchase. Written authorization must be obtained from Micro-Tech Endoscopy prior to return.

PRODUCT STORAGE

Micro-Tech Endoscopy products should be stored in a cool, dry, clean, well-ventilated place. Do not expose the packaging to organic solvents, ionizing radiation or ultraviolet radiation.

DESIGN CHANGES

Micro-Tech Endoscopy reserves the right to discontinue or change products without notice. If any item has been discontinued, alternatives will be offered.

UPN DESCRIPTION

Now products may be ordered using the UPN (Order Number) as well as the traditional Reference Number (Catalogue Number).

QUALITY POLICY

Micro-Tech Endoscopy is committed to supplying a comprehensive portfolio of quality, cost-effective gastrointestinal devices. We are committed to maintaining a quality system that provides for safe and effective products and services that meet the needs of our patients, customers and stakeholders.

MICRO-TECH®
ENDOSCOPY

US.MTMED.COM

©2026 Micro-Tech Endoscopy USA Inc. All rights reserved.
All other trademarks are the property of their respective owners.
MPW10000H Rev 6